

Factor IX19. GlycoPEGylation of Factor IX produced in CHO cells

This example sets forth the preparation of asialoFactor IX and its sialylation with CMP-sialic acid-PEG.

- 5       **Desialylation of rFactor IX.** A recombinant form of Coagulation Factor IX (rFactor IX ) was made in CHO cells. 6000 IU of rFactor IX were dissolved in a total of 12 mL USP H<sub>2</sub>O. This solution was transferred to a Centricon Plus 20, PL-10 centrifugal filter with another 6 mL USP H<sub>2</sub>O. The solution was concentrated to 2 mL and then diluted with 15 mL 50 mM Tris-HCl pH 7.4, 0.15 M NaCl, 5 mM CaCl<sub>2</sub>, 0.05% NaN<sub>3</sub> and then reconcentrated.
- 10      The dilution/concentration was repeated 4 times to effectively change the buffer to a final volume of 3.0 mL. Of this solution, 2.9 mL (about 29 mg of rFactor IX) was transferred to a small plastic tube and to it was added 530 mU  $\alpha$ 2-3,6,8-Neuraminidase– agarose conjugate (*Vibrio cholerae*, Calbiochem, 450  $\mu$ L). The reaction mixture was rotated gently for 26.5 hours at 32 °C. The mixture was centrifuged 2 minutes at 10,000 rpm and the supernatant
- 15      was collected. The agarose beads (containing neuraminidase) were washed 6 times with 0.5 mL 50 mM Tris-HCl pH 7.12, 1 M NaCl, 0.05% NaN<sub>3</sub>. The pooled washings and supernatants were centrifuged again for 2 minutes at 10,000 rpm to remove any residual agarose resin. The pooled, desialylated protein solution was diluted to 19 mL with the same buffer and concentrated down to ~ 2 mL in a Centricon Plus 20 PL-10 centrifugal filter. The
- 20      solution was twice diluted with 15 mL of 50 mM Tris-HCl pH 7.4, 0.15 M NaCl, 0.05% NaN<sub>3</sub> and reconcentrated to 2 mL. The final desialylated rFactor IX solution was diluted to 3 mL final volume (~10 mg/mL) with the Tris Buffer. Native and desialylated rFactor IX samples were analyzed by IEF-Electrophoresis. Isoelectric Focusing Gels (pH 3-7) were run using 1.5  $\mu$ L (15  $\mu$ g) samples first diluted with 10  $\mu$ L Tris buffer and mixed with 12  $\mu$ L
- 25      sample loading buffer. Gels were loaded, run and fixed using standard procedures. Gels were stained with Colloidal Blue Stain (Figure 154), showing a band for desialylated Factor IX.

- Preparation of PEG (1 kDa and 10 kDa)-SA-Factor IX.** Desialylated rFactor-IX (29 mg, 3 mL) was divided into two 1.5 mL (14.5 mg) samples in two 15 mL centrifuge
- 30      tubes. Each solution was diluted with 12.67 mL 50 mM Tris-HCl pH 7.4, 0.15 M NaCl, 0.05% NaN<sub>3</sub> and either CMP-SA-PEG-1k or 10k (7.25  $\mu$ mol) was added. The tubes were

inverted gently to mix and 2.9 U ST3Gal3 (326  $\mu$ L) was added (total volume 14.5 mL). The tubes were inverted again and rotated gently for 65 hours at 32 °C. The reactions were stopped by freezing at -20 °C. 10  $\mu$ g samples of the reactions were analyzed by SDS-PAGE. The PEGylated proteins were purified on a Toso Haas Biosep G3000SW (21.5 x 30 cm, 13  $\mu$ m) HPLC column with Dulbecco's Phosphate Buffered Saline, pH 7.1 (Gibco), 6 mL/min. The reaction and purification were monitored using SDS Page and IEF gels. Novex Tris-Glycine 4-20% 1 mm gels were loaded with 10  $\mu$ L (10  $\mu$ g) of samples after dilution with 2  $\mu$ L of 50 mM Tris-HCl, pH 7.4, 150 mM NaCl, 0.05% NaN<sub>3</sub> buffer and mixing with 12  $\mu$ L sample loading buffer and 1  $\mu$ L 0.5 M DTT and heated for 6 minutes at 85 °C. Gels were stained with Colloidal Blue Stain (Figure 155) showing a band for PEG (1 kDa and 10 kDa)-SA-Factor IX.

## 20. Direct Sialyl-GlycoPEGylation of Factor IX

This example sets forth the preparation of sialyl-PEGylation of Factor IX without prior sialidase treatment.

**Sialyl-PEGylation of Factor-IX with CMP-SA-PEG-(10 kDa).** Factor IX (1100 IU), which was expressed in CHO cells and was fully sialylated, was dissolved in 5 mL of 20 mM histidine, 520 mM glycine, 2% sucrose, 0.05% NaN<sub>3</sub> and 0.01% polysorbate 80, pH 5.0. The CMP-SA-PEG-(10 kDa) (27 mg, 2.5  $\mu$ mol) was then dissolved in the solution and 1 U of ST3Gal3 was added. The reaction was complete after gently mixing for 28 hours at 32°C. The reaction was analyzed by SDS-PAGE as described by Invitrogen. The product protein was purified on an Amersham Superdex 200 (10 x 300 mm, 13  $\mu$ m) HPLC column with phosphate buffered saline, pH 7.0 (PBS), 1 mL/min.  $R_t$  = 9.5 min.

**Sialyl-PEGylation of Factor-IX with CMP-SA-PEG-(20 kDa).** Factor IX (1100 IU), which was expressed in CHO cells and was fully sialylated, was dissolved in 5 mL of 20 mM histidine, 520 mM glycine, 2% sucrose, 0.05% NaN<sub>3</sub> and 0.01% polysorbate 80, pH 5.0. The CMP-SA-PEG-(20 kDa) (50 mg, 2.3  $\mu$ mol) was then dissolved in the solution and CST-II was added. The reaction mixture was complete after gently mixing for 42 hours at 32°C. The reaction was analyzed by SDS-PAGE as described by Invitrogen.

The product protein was purified on an Amersham Superdex 200 (10 x 300 mm, 13  $\mu$ m) HPLC column with phosphate buffered saline, pH 7.0 (Fisher), 1 mL/min.  $R_t$  = 8.6 min.

## 21. Sialic Acid Capping of GlycoPEGylated Factor IX

5 This examples sets forth the procedure for sialic acid capping of sialyl-glycoPEGylated peptides. Here, Factor-IX is the exemplary peptide.

**Sialic acid capping of N-linked and O-linked Glycans of Factor-IX-SA-PEG (10 kDa).** Purified r-Factor-IX-PEG (10 kDa) (2.4 mg) was concentrated in a Centricon<sup>®</sup> Plus 20 PL-10 (Millipore Corp., Bedford, MA) centrifugal filter and the buffer was changed to 50  
10 mM Tris-HCl pH 7.2, 0.15 M NaCl, 0.05% NaN<sub>3</sub> to a final volume of 1.85 mL. The protein solution was diluted with 372  $\mu$ L of the same Tris buffer and 7.4 mg CMP-SA (12  $\mu$ mol) was added as a solid. The solution was inverted gently to mix and 0.1 U ST3Gal1 and 0.1 U ST3Gal3 were added. The reaction mixture was rotated gently for 42 hours at 32 °C.

A 10  $\mu$ g sample of the reaction was analyzed by SDS-PAGE. Novex Tris-Glycine 4-  
15 12% 1 mm gels were performed and stained using Colloidal Blue as described by Invitrogen. Briefly, samples, 10  $\mu$ L (10  $\mu$ g), were mixed with 12  $\mu$ L sample loading buffer and 1  $\mu$ L 0.5 M DTT and heated for 6 minutes at 85 °C (Figure 156, lane 4).

## Factor VIIa

### 20 22. GlycoPEGylation of Recombinant Factor VIIa produced in BHK cells

This example sets forth the PEGylation of recombinant Factor VIIa made in BHK cells.

**Preparation of Asialo-Factor VIIa.** Recombinant Factor VIIa was produced in BHK cells (baby hamster kidney cells). Factor VIIa (14.2 mg) was dissolved at 1 mg/ml in  
25 buffer solution (pH 7.4, 0.05 M Tris, 0.15 M NaCl, 0.001 M CaCl<sub>2</sub>, 0.05% NaN<sub>3</sub>) and was incubated with 300 mU/mL sialidase (*Vibrio cholera*)-agarose conjugate for 3 days at 32 °C. To monitor the reaction a small aliquot of the reaction was diluted with the appropriate buffer and an IEF gel performed according to Invitrogen procedures (Figure 157). The mixture was centrifuged at 3,500 rpm and the supernatant was collected. The resin was washed three  
30 times (3 $\times$ 2 mL) with the above buffer solution ( pH 7.4, 0.05 M Tris, 0.15 M NaCl, 0.05% NaN<sub>3</sub>) and the combined washes were concentrated in a Centricon-Plus-20. The remaining

solution was buffer exchanged with 0.05 M Tris (pH 7.4), 0.15 M NaCl, 0.05% NaN<sub>3</sub> to a final volume of 14.4 mL.

**Preparation of Factor VIIa-SA-PEG (1 kDa and 10 kDa).** The desialylation of Factor VIIa solution was split into two equal 7.2 ml samples. To each sample was added either CMP-SA-5-PEG(1 kDa) (7.4 mg) or CMP-SA-5-PEG(10 kDa) (7.4 mg). ST3Gal3 (1.58U) was added to both tubes and the reaction mixtures were incubated at 32°C for 96 hrs. The reaction was monitored by SDS-PAGE gel using reagents and conditions described by Invitrogen. When the reaction was complete, the reaction mixture was purified using a Toso Haas TSK-Gel-3000 preparative column using PBS buffer (pH 7.1) and collecting fractions based on UV absorption. The combined fractions containing the product were concentrated at 4°C in Centricon-Plus-20 centrifugal filters (Millipore, Bedford, MA) and the concentrated solution reformulated to yield 1.97 mg (bicinchoninic acid protein assay, BCA assay, Sigma-Aldrich, St. Louis MO) of Factor VIIa-PEG. The product of the reaction was analyzed using SDS-PAGE and IEF analysis according to the procedures and reagents supplied by Invitrogen. Samples were dialyzed against water and analyzed by MALDI-TOF. Figure 158 shows the MALDI results for native Factor VIIa. Figure 159 contains the MALDI results for Factor VIIa PEGylated with 1 kDa PEG where peak of Factor VIIa PEGylated with 1KDa PEG is evident. Figure 160 contains the MALDI results for Factor VIIa PEGylated with 10 kDa PEG where a peak for Factor VIIa PEGylated with 10 kDa PEG is evident. Figure 161 depicts the SDS-PAGE analysis of all of the reaction products, where a band for Factor VIIa-SA-PEG (10 kDa) is evident.

### Follicle Stimulating Hormone (FSH)

#### 23. GlycoPEGylation of human pituitary-derived FSH

This example illustrates the assembly of a conjugate of the invention. Follicle Stimulating Hormone (FSH) is desialylated and then conjugated with CMP-(sialic acid)-PEG.

**Desialylation of Follicle Stimulating Hormone.** Follicle Stimulating Hormone (FSH) (Human Pituitary, Calbiochem Cat No. 869001), 1 mg, was dissolved in 500 µL 50 mM Tris-HCl pH 7.4, 0.15 M NaCl, 5 mM CaCl<sub>2</sub>. This solution, 375 µL, was transferred to a small plastic tube and to it was added 263 mU Neuraminidase II (*Vibrio cholerae*). The reaction mixture was shaken gently for 15 hours at 32 °C. The reaction mixture was added to

N-(*p*-aminophenyl)oxamic acid-agarose conjugate, 600  $\mu$ L, pre-equilibrated with 50 mM Tris-HCl pH 7.4, 150 mM NaCl and 0.05% NaN<sub>3</sub> and gently rotated 6.5 hours at 4 °C. The suspension was centrifuged for 2 minutes at 14,000 rpm and the supernatant was collected. The beads were washed 5 times with 0.5 mL of the buffer and all supernatants were pooled.

5 The enzyme solution was dialyzed (7000 MWCO) for 15 hours at 4 °C with 2 L of a solution containing 50 mM Tris -HCl pH 7.4, 1 M NaCl, 0.05% NaN<sub>3</sub>, and then twice for 4 hours at 4 °C into 50 mM Tris -HCl pH 7.4, 1 M NaCl, 0.05% NaN<sub>3</sub>. The solution was concentrated to 2  $\mu$ g/ $\mu$ L by Speed Vac and stored at -20 °C. Reaction samples were analyzed by IEF gels (pH 3-7) (Invitrogen) (Figure 162).

10 **Preparation of human pituitary-derived SA-FSH and PEG-SA-Follicle Stimulating Hormone.** Desialylated FSH (100  $\mu$ g, 50  $\mu$ L) and CMP-sialic acid or CMP-SA-PEG (1 kDa or 10 kDa) (0.05  $\mu$ mol) were dissolved in 13.5  $\mu$ L H<sub>2</sub>O (adjusted to pH 8 with NaOH) in 0.5 mL plastic tubes. The tubes were vortexed briefly and 40 mU ST3Gal3 (36.5  $\mu$ L) was added (total volume 100  $\mu$ L). The tubes were vortexed again and shaken gently for

15 24 hours at 32 °C. The reactions were stopped by freezing at -80 °C. Reaction samples of 15  $\mu$ g were analyzed by SDS-PAGE (Figure 163), IEF gels (Figure 164) and MALDI-TOF. Native FSH was also analyzed by SDS-PAGE (Figure 165)

**Analysis of SDS PAGE and IEF Gels of Reaction Products.** Novex Tris-Glycine 8-16% 1 mm gels for SDS PAGE analysis were purchased from Invitrogen. 7.5  $\mu$ L (15  $\mu$ g)

20 of FSH reaction samples were diluted with 5  $\mu$ L of 50 mM Tris-HCl, pH 7.4, 150 mM NaCl, 0.05% NaN<sub>3</sub> buffer, mixed with 15  $\mu$ L sample loading buffer and 1  $\mu$ L 9 M  $\mu$ -mercaptoethanol and heated for 6 minutes at 85 °C. Gels were run as directed by Invitrogen and stained with Colloidal Blue Stain (Invitrogen).

FSH samples (15  $\mu$ g) were diluted with 5  $\mu$ L Tris buffer and mixed with 15  $\mu$ L

25 sample loading buffer (Figure 162). The samples were then applied to Isoelectric Focusing Gels (pH 3-7) (Invitrogen) (Figure 165). Gels were run and fixed as directed by Invitrogen and then stained with Colloidal Blue Stain.

24. GlycoPEGylation of recombinant FSH produced recombinantly in CHO cells

This example illustrates the assembly of a conjugate of the invention. Desialylated FSH was conjugated with CMP-(sialic acid)-PEG.

5       **Preparation of recombinant Asialo-Follicle Stimulation Hormone.** Recombinant Follicle Stimulation Hormone (rFSH) produced from CHO was used in these studies. The 7,500 IU of rFSH was dissolved in 8 mL of water. The FSH solution was dialyzed in 50 mM Tris-HCl pH 7.4, 0.15 M NaCl, 5 mM CaCl<sub>2</sub> and concentrated to 500 µL in a Centricon Plus 20 centrifugal filter. A portion of this solution (400 µL) (~ 0.8 mg FSH) was transferred to a  
10   small plastic tube and to it was added 275 mU Neuraminidase II (*Vibrio cholerae*). The reaction mixture was mixed for 16 hours at 32 °C. The reaction mixture was added to prewashed N-(*p*-aminophenyl)oxamic acid-agarose conjugate (800 µL) and gently rotated for 24 hours at 4 °C. The mixture was centrifuged at 10,000 rpm and the supernatant was collected. The beads were washed 3 times with 0.6 mL Tris-EDTA buffer, once with 0.4 mL  
15   Tris-EDTA buffer and once with 0.2 mL of the Tris-EDTA buffer and all supernatants were pooled. The supernatant was dialyzed at 4 °C against 2 L of 50 mM Tris -HCl pH 7.4, 1 M NaCl, 0.05% NaN<sub>3</sub> and then twice more against 50 mM Tris -HCl pH 7.4, 1 M NaCl, 0.05% NaN<sub>3</sub>. The dialyzed solution was then concentrated to 420 µL in a Centricon Plus 20 centrifugal filter and stored at -20 °C.

20       Native and desialylated rFSH samples were analyzed by SDS-PAGE and IEF (Figure 166). Novex Tris-Glycine 8-16% 1 mm gels were purchased from Invitrogen. Samples (7.5 µL, 15 µg) samples were diluted with 5 µL of 50 mM Tris-HCl, pH 7.4, 150 mM NaCl, 0.05% NaN<sub>3</sub> buffer, mixed with 15 µL sample loading buffer and 1 µL 9 M β-mercaptoethanol and heated for 6 minutes at 85 °C. Gels were run as directed by Invitrogen  
25   and stained with Colloidal Blue Stain (Invitrogen). Isoelectric Focusing Gels (pH 3-7) were purchased from Invitrogen. Samples (7.5 µL, 15 µg) were diluted with 5 µL Tris buffer and mixed with 15 µL sample loading buffer. Gels were loaded, run and fixed as directed by Invitrogen. Gels were stained with Colloidal Blue Stain. Samples of native and desialylated FSH were also dialyzed against water and analyzed by MALDI-TOF.

30       **Sialyl-PEGylation of recombinant Follicle Stimulation Hormone.** Desialylated FSH (100 µg, 54 µL) and CMP-SA-PEG (1 kDa or 10 kDa) (0.05 µmol) were dissolved in 28

μL 50 mM Tris-HCl, 0.15 M NaCl, 0.05% NaN<sub>3</sub>, pH 7.2 in 0.5 mL plastic tubes. The tubes were vortexed briefly and 20 mU of ST3Gal3 was added (total volume 100 μL). The tubes were vortexed again, mixed gently for 24 hours at 32 °C and the reactions stopped by freezing at -80 °C. Samples of this reaction were analyzed as described above by SDS-PAGE gels (Figure 167), IEF gels (Figure 168) and MALDI-TOF MS.

MALDI was also performed on the PEGylated rFSH. During ionization, SA-PEG is eliminated from the N-glycan structure of the glycoprotein. Native FSH gave a peak at 13928; AS-rFSH (13282); resialylated r-FSH (13332); PEG1000-rFSH (13515; 14960 (1); 16455 (2); 17796 (3); 19321 (4)); and PEG 10000 (23560 (1); 34790 (2); 45670 (3); and 56760 (4)).

#### 25. Pharmacokinetic Study of GlycoPEGylated FSH

This example sets forth the *in vivo* testing of the pharmacokinetic properties glycoPEGylated Follicle Stimulating Hormone (FSH) prepared according to the methods of the invention as compared to non-PEGylated FSH.

FSH, FSH-SA-PEG (1 kDa) and FSH-SA-PEG (10 kDa) were radioiodinated using standard conditions (Amersham Biosciences, Arlington Heights, IL) and formulated in phosphate buffered saline containing 0.1% BSA. After dilution in phosphate buffer to the appropriate concentration, each of the test FSH proteins (0.4 μg, each) was injected intravenously into female Sprague Dawley rats (250-300 g body weight) and blood drawn at time points from 0 to 80 hours. Radioactivity in blood samples was analyzed using a gamma counter and the pharmacokinetics analyzed using standard methods (Figure 169). FSH was cleared from the blood much more quickly than FSH-PEG(1 kDa), which in turn was clear somewhat more quickly than FSH-PEG(10 kDa).

#### 26. Sertoli Cell Bioassay for *In Vitro* Activity of GlycoPEGylated FSH

This example sets forth a bioassay for follicle stimulating hormone (FSH) activity based on cultured Sertoli cells. This assay is useful to determine the bioactivity of FSH after glycan remodeling, including glycoconjugation.

This bioassay is based on the dose-response relationship that exists between the amount of estradiol produced when FSH, but not lutenizing hormone (LH), is added to

cultured Sertoli cells obtained from immature old rats. Exogenous testosterone is converted to 17 $\beta$ -estradiol in the presence of FSH.

Seven to 10 days old Sprague-Dawley rats were used to obtain Sertoli cells. After sacrifice, testes were decapsulated and tissue was dispersed by incubation in collagenase (1 mg/ml), trypsin (1mg/ml), hyaluronidase (1 mg/ml) and DNases (5  $\mu$ g/ml) for 5 to 10 min. The tubule fragments settled to the bottom of the flask and were washed in PBS (1x). The tubule fragments were reincubated for 20 min with a media containing the same enzymes: collagenase (1 mg/ml), trypsin (1mg/ml), hyaluronidase (1 mg/ml) and DNases (5  $\mu$ g/ml).

The tubule fragments were homogenized and plated into a 24 well plate in a serum free media. 5 x 10<sup>5</sup> cells were dispersed per well. After 48h incubation at 37° C and 5% CO<sub>2</sub>, fresh media was added to the cells. Composition of the serum free media: DMEM (1 vol), Ham's F10 nutrient mixture (1 vol), insulin 1  $\mu$ g/ml, Transferrin 5  $\mu$ g/ml, EGF 10 ng/ml, T4 20 pg/ml, Hydrocortisone 10<sup>-8</sup> M, Retinoic acid 10<sup>-6</sup> M.

The stimulation experiment consists of a 24 hour incubation with standard FSH or samples at 37°C and 5% CO<sub>2</sub>. The mean intra-assay coefficient of variation is 9% and the mean inter-assay coefficient of variation is 11%.

The 17B-estradiol Elisa Kit DE2000 (R&D Systems, Minneapolis, MN) was used to quantify the level of estradiol after incubation with FSH, FSH-SA-PEG (1 kDa) and FSH-SA-PEG (10 kDa).

The procedure was as follows: 100  $\mu$ l of Estradiol Standard (provided with kit and prepared as per instructions with kit) or sample was pipetted into wells of 17B-estradiol Elisa plate(s); 50  $\mu$ l of 17B-estradiol Conjugate (provided with kit, prepared as per instructions with kit) was added to each well; 50  $\mu$ l of 17B-estradiol antibody solution (provided with kit and prepared as per instructions with kit) was added to each well; plates were incubated for 2 hour at room temperature at 200 rpm; the liquid was aspirated from each well; the wells were washed 4 times using the washing solution; all the liquid was removed from the wells; 200  $\mu$ l of pNPP Substrate (provided with kit and prepared as per instructions with kit) was added to all wells and incubated for 45 min; 50  $\mu$ l of Stop solution (provided with kit and prepared as per instructions with kit) was added and the plates were read it at 405 nm (Figure 170).

While FSH-PEG(10 kDa) exhibited a modest stimulation of Sertoli cells, at 1  $\mu$ g/ml, FSH-PEG(1 kDa) stimulated Sertoli cells up to 50% more than unPEGylated FSH.



### 27. Steelman-Pohley Bioassay of *In Vivo* Activity of GlycoPEGylated FSH

In this example, the Steelman-Pohley bioassay (Steelman and Pohley, 1953, Endocrinology 53:604-615) was used to determine the *in vivo* activity of glycoPEGylated FSH. The Steelman-Pohley assay uses the change in ovary weight of a rat to measure the *in vivo* activity of FSH that is coinjected with human chorionic gonadotropin.

The Steelman-Pohley bioassay was performed according to the protocol described in Christin-Maitre et al. (2000, Methods 21:51-57). Seventy female Sprague-Dawley Rats (Charles River Laboratories, Wilmington, MA), aged 21 to 22 days, were housed in the testing facility for at least 5 days before the beginning the assay procedure. Throughout the procedure, the animal room was climate controlled at 18 to 26°C, 30 to 70% relative humidity, and 12 hr. artificial light/12 hr. dark. All animals were fed Certified Rodent Chow (Harlan Teklad, Madison WI) or the equivalent, and water, both *ad libitum*. Animal procedures were performed at Calvert Preclinical Services, Inc. (Olyphant, PA).

Recombinant FSH was expressed in CHO cells, purified by standard techniques and glycoPEGylated with PEG (1 kDa). The rats were divided into seven test groups, with ten animals per group. On days -1 and 0, animals of all groups were subcutaneously injected with 20 I.U. of human chorionic gonadotropin (HCG) in 0.5 ml of 0.9 % NaCl. On days 1, 2 and 3, the control animals were subcutaneously injected with a dose of 0.5 ml containing 20 I.U. HCG in 0.9% NaCl, while in the other groups, the HCG dose was augmented with either rFSH or rFSH-SA-PEG (1 kDa) at either 0.14 µg, 0.4 µg or 1.2 µg per dose. On day 4, the animals were euthanized by CO<sub>2</sub> inhalation. The ovaries were removed, trimmed and weighted. The average ovary weight was determined for each group.

Figure 171 presents the average ovary weight of the test groups on day 4. The groups receiving HCG alone (control) or the low dose (0.14 µg) of either rFSH or rFSH-SA-PEG (1 kDa) had ovary weights that were roughly equivalent. The groups receiving the medium (0.4 µg) or high (1.2 µg) doses of rFSH or rFSH-SA-PEG (1 kDa) had ovary weights roughly twice that of the control group. At the medium dose (0.4 µg), the glycoPEGylated rFSH had roughly the same *in vivo* activity (as determined by ovary weight) as the unPEGylated rFSH.

At the high dose (1.2 µg), the glycoPEGylated rFSH had somewhat higher *in vivo* activity than the unPEGylated rFSH.

### G-CSF

#### 5                   28. GlycoPEGylation of G-CSF produced in CHO cells

**Preparation of Asialo-Granulocyte-Colony Stimulation Factor (G-CSF).** G-CSF produced in CHO cells is dissolved at 2.5 mg/mL in 50 mM Tris 50 mM Tris-HCl pH 7.4, 0.15 M NaCl, 5 mM CaCl<sub>2</sub> and concentrated to 500 µL in a Centricon Plus 20 centrifugal filter. The solution is incubated with 300 mU/mL Neuraminidase II (*Vibrio cholerae*) for 16  
10   hours at 32 °C. To monitor the reaction a small aliquot of the reaction is diluted with the appropriate buffer and a IEF gel performed. The reaction mixture is then added to prewashed N-(*p*-aminophenyl)oxamic acid-agarose conjugate (800 µL/mL reaction volume) and the washed beads gently rotated for 24 hours at 4 °C. The mixture is centrifuged at 10,000 rpm and the supernatant was collected. The beads are washed 3 times with Tris-EDTA buffer,  
15   once with 0.4 mL Tris-EDTA buffer and once with 0.2 mL of the Tris-EDTA buffer and all supernatants are pooled. The supernatant is dialyzed at 4 °C against 50 mM Tris –HCl pH 7.4, 1 M NaCl, 0.05% NaN<sub>3</sub> and then twice more against 50 mM Tris –HCl pH 7.4, 1 M NaCl, 0.05% NaN<sub>3</sub>. The dialyzed solution is then concentrated using a Centricon Plus 20 centrifugal filter and stored at –20 °C. The conditions for the IEF gel were run according to  
20   the procedures and reagents provided by Invitrogen. Samples of native and desialylated G-CSF are dialyzed against water and analyzed by MALDI-TOF MS.

**Preparation of G-CSF-(alpha2,3)-Sialyl-PEG.** Desialylated G-CSF was dissolved at 2.5 mg/mL in 50 mM Tris-HCl, 0.15 M NaCl, 0.05% NaN<sub>3</sub>, pH 7.2. The solution is incubated with 1 mM CMP-sialic acid-PEG and 0.1 U/mL of ST3Gal1 at 32°C for 2 days.  
25   To monitor the incorporation of sialic acid-PEG, a small aliquot of the reaction had CMP-SA-PEG-fluorescent ligand added; the label incorporated into the peptide is separated from the free label by gel filtration on a Toso Haas G3000SW analytical column using PBS buffer (pH 7.1). The fluorescent label incorporation into the peptide is quantitated using an in-line fluorescent detector. After 2 days, the reaction mixture is purified using a Toso Haas  
30   G3000SW preparative column using PBS buffer (pH 7.1) and collecting fractions based on UV absorption. The product of the reaction is analyzed using SDS-PAGE and IEF analysis

according to the procedures and reagents supplied by Invitrogen. Samples of native and PEGylated G-CSF are dialyzed against water and analyzed by MALDI-TOF MS.

**Preparation of G-CSF-(alpha2,8)-Sialyl-PEG.** G-CSF produced in CHO cells, which contains an alpha2,3-sialylated O-linked glycan, is dissolved at 2.5 mg/mL in 50 mM Tris-HCl, 0.15 M NaCl, 0.05% NaN<sub>3</sub>, pH 7.2. The solution is incubated with 1 mM CMP-sialic acid-PEG and 0.1 U/mL of CST-II at 32°C for 2 days. To monitor the incorporation of sialic acid-PEG, a small aliquot of the reaction has CMP-SA-PEG-fluorescent ligand added; the label incorporated into the peptide is separated from the free label by gel filtration on a Toso Haas G3000SW analytical column using PBS buffer (pH 7.1). The fluorescent label incorporation into the peptide is quantitated using an in-line fluorescent detector. After 2 days, the reaction mixture is purified using a Toso Haas G3000SW preparative column using PBS buffer (pH 7.1) and collecting fractions based on UV absorption. The product of the reaction is analyzed using SDS-PAGE and IEF analysis according to the procedures and reagents supplied by Invitrogen. Samples of native and PEGylated G-CSF are dialyzed against water and analyzed by MALDI-TOF MS.

**Preparation of G-CSF-(alpha2,6)-Sialyl-PEG.** G-CSF, containing only O-linked GalNAc, is dissolved at 2.5 mg/mL in 50 mM Tris-HCl, 0.15 M NaCl, 0.05% NaN<sub>3</sub>, pH 7.2. The solution is incubated with 1 mM CMP-sialic acid-PEG and 0.1 U/mL of ST6GalNAcI or II at 32°C for 2 days. To monitor the incorporation of sialic acid-PEG, a small aliquot of the reaction has CMP-SA-PEG-fluorescent ligand added; the label incorporated into the peptide is separated from the free label by gel filtration on a Toso Haas G3000SW analytical column using PBS buffer (pH 7.1). The fluorescent label incorporation into the peptide is quantitated using an in-line fluorescent detector. After 2 days, the reaction mixture is purified using a Toso Haas G3000SW preparative column using PBS buffer (pH 7.1) and collecting fractions based on UV absorption. The product of the reaction is analyzed using SDS-PAGE and IEF analysis according to the procedures and reagents supplied by Invitrogen. Samples of native and PEGylated G-CSF are dialyzed against water and analyzed by MALDI-TOF MS.

G-CSF produced in CHO cells was treated with Arthrobacter sialidase and was then purified by size exclusion on Superdex75 and was treated with ST3Gal1 or ST3 Gal2 and then with CMP-SA-PEG 20Kda. The resulting molecule was purified by ion exchange and

gel filtration and analysis by SDS PAGE demonstrated that the PEGylation was complete. This is the first demonstration of glycoPEGylation of an O-linked glycan.

### Glucocerebrosidase

#### 5                    29. Glucocerebrosidase-mannose-6-phosphate produced in CHO cells

This example sets forth the procedure to glycoconjugate mannose-6-phosphate to a peptide produced in CHO cells such as glucocerebrosidase.

**Preparation of asialo-glucoceramide.** Glucocerebrosidase produced in CHO cells is dissolved at 2.5 mg/mL in 50 mM Tris 50 mM Tris-HCl pH 7.4, 0.15 M NaCl, and is  
10 incubated with 300 mU/mL sialidase-agarose conjugate for 16 hours at 32 °C. To monitor the reaction a small aliquot of the reaction is diluted with the appropriate buffer and a IEF gel and SDS-PAGE performed according to Invitrogen procedures. The mixture is centrifuged at 10,000 rpm and the supernatant is collected. The beads are washed 3 times with Tris-EDTA buffer, once with 0.4 mL Tris-EDTA buffer, and once with 0.2 mL of the Tris-EDTA buffer.  
15 All supernatants are pooled. The supernatant is dialyzed at 4 °C against 50 mM Tris-HCl pH 7.4, 1 M NaCl, 0.05% NaN<sub>3</sub> and then twice more against 50 mM Tris-HCl pH 7.4, 1 M NaCl, 0.05% NaN<sub>3</sub>. The dialyzed solution is then concentrated using a Centricon Plus 20 centrifugal filter. The product of the reaction is analyzed using SDS-PAGE and IEF analysis according to the procedures and reagents supplied by Invitrogen. Samples are dialyzed  
20 against water and analyzed by MALDI-TOF MS.

#### **Preparation of Glucocerebrosidase-SA-linker-Mannose-6-phosphate (procedure**

**1).** Asialo-glucocerebrosidase from above is dissolved at 2.5 mg/mL in 50 mM Tris-HCl, 0.15 M NaCl, 0.05% NaN<sub>3</sub>, pH 7.2. The solution is incubated with 1 mM CMP-sialic acid-linker-Man-6-phosphate and 0.1 U/mL of ST3Gal3 at 32°C for 2 days. To monitor the  
25 incorporation of sialic acid-linker-Man-6-phosphate, a small aliquot of the reaction had CMP-SA-PEG-fluorescent ligand added; the label incorporated into the peptide is separated from the free label by gel filtration on a Toso Haas TSK-Gel-3000 analytical column using PBS buffer (pH 7.1). The fluorescent label incorporation into the peptide is quantitated using an in-line fluorescent detector. When the reaction is complete, the reaction mixture is purified  
30 using a Toso Haas TSK-Gel-3000 preparative column using PBS buffer (pH 7.1) and collecting fractions based on UV absorption. The product of the reaction is analyzed using

SDS-PAGE and IEF analysis according to the procedures and reagents supplied by Invitrogen. Samples are dialyzed against water and analyzed by MALDI-TOF MS.

**Preparation of Glucocerebrosidase-SA-linker-Mannose-6-phosphate (procedure 2).** Glucocerebrosidase, produced in CHO but incompletely sialylated, is dissolved at 2.5 mg/mL in 50 mM Tris-HCl, 0.15 M NaCl, 0.05% NaN<sub>3</sub>, pH 7.2. The solution is incubated with 1 mM CMP-sialic acid-linker-Man-6-phosphate and 0.1 U/mL of ST3Gal3 at 32°C for 2 days. To monitor the incorporation of sialic acid-linker-Man-6-phosphate, a small aliquot of the reaction had CMP-SA-PEG-fluorescent ligand added; the label incorporated into the peptide is separated from the free label by gel filtration on a Toso Haas TSK-Gel-3000 analytical column using PBS buffer (pH 7.1). The fluorescent label incorporation into the peptide is quantitated using an in-line fluorescent detector. When the reaction is complete, the reaction mixture is purified using a Toso Haas TSK-Gel-3000 preparative column using PBS buffer (pH 7.1) and collecting fractions based on UV absorption. The product of the reaction is analyzed using SDS-PAGE and IEF analysis according to the procedures and reagents supplied by Invitrogen. Samples are dialyzed against water and analyzed by MALDI-TOF MS.

### 30. Glucocerebrosidase-transferrin

This example sets forth the procedures for the glycoconjugation of proteins, and in particular, transferrin is glycoconjugated to glucocerebrosidase. The GlcNAc-ASN structures are created on glucoceraminidase, and Transferrin-SA-Linker-Gal-UDP is conjugated to GNDF GlcNAc-ASN structures using galactosyltransferase.

**Preparation of GlcNAc-glucocerebrosidase (Cerezyme™).** Cerezyme™ (glucocerebrosidase) produced in CHO cells is dissolved at 2.5 mg/mL in 50 mM Tris 50 mM Tris-HCl pH 7.4, 0.15 M NaCl, and is incubated with 300 mU/mL Endo-H-agarose conjugate for 16 hours at 32 °C. To monitor the reaction a small aliquot of the reaction is diluted with the appropriate buffer and a IEF gel and SDS-PAGE performed according to Invitrogen procedures. The mixture is centrifuged at 10,000 rpm and the supernatant is collected. The beads are washed 3 times with Tris-EDTA buffer, once with 0.4 mL Tris-EDTA buffer and once with 0.2 mL of the Tris-EDTA buffer and all supernatants are pooled. The supernatant is dialyzed at 4 °C against 50 mM Tris -HCl pH 7.4, 1 M NaCl, 0.05% NaN<sub>3</sub> and then twice

more against 50 mM Tris -HCl pH 7.4, 1 M NaCl, 0.05% NaN<sub>3</sub>. The dialyzed solution is then concentrated using a Centricon Plus 20 centrifugal filter. The product of the reaction is analyzed using SDS-PAGE and IEF analysis according to the procedures and reagents supplied by Invitrogen. Samples are dialyzed against water and analyzed by MALDI-TOF MS.

**Preparation of Transferrin-SA-Linker-Gal-glucocerebrosidase.** Transferrin-SA-Linker-Gal-UDP from above is dissolved at 2.5 mg/mL in 50 mM Tris-HCl, 0.15 M NaCl, 5 mM MnCl<sub>2</sub>, 0.05% NaN<sub>3</sub>, pH 7.2. The solution is incubated with 2.5 mg/mL GlcNAc-glucocerebrosidase and 0.1 U/mL of galactosyltransferase at 32°C for 2 days. To monitor the incorporation of glucocerebrosidase, the peptide is separated by gel filtration on a Toso Haas G3000SW analytical column using PBS buffer (pH 7.1) and the product detected by UV absorption. The reaction mixture is then purified using a Toso Haas G3000SW preparative column using PBS buffer (pH 7.1) collecting fractions based on UV absorption. The product of the reaction is analyzed using SDS-PAGE and IEF analysis according to the procedures and reagents supplied by Invitrogen. Samples are dialyzed against water and analyzed by MALDI-TOF MS.

#### GM-CSF

##### 31. Generation and PEGylation of GlcNAc-ASN Structures: GM-CSF produced in *Saccharomyces*

This example sets forth the preparation of Tissue-type Activator with PEGylated GlcNAc-Asn structures.

Recombinant GM-CSF expressed in yeast is expected to contain 2 N-linked and 2 O-linked glycans. The N-linked glycans should be of the branched mannan type. This recombinant glycoprotein is treated with an endoglycosidase from the group consisting of endoglycosidase H, endoglycosidase-F1, endoglycosidase-F2, endoglycosidase-F3, endoglycosidase-M either alone or in combination with mannosidases I, II and III to generate GlcNAc nubs on the asparagine (Asn) residues on the peptide/protein backbone.

The GlcNAc-Asn structures on the peptide/protein backbone is then be modified with galactose or galactose-PEG using UDP-galactose or UDP-galactose-6-PEG, respectively, and a galactosyltransferase such as GalT1. In one case the galactose-PEG is the terminal residue.

In the second case the galactose is further modified with SA-PEG using a CMP-SA-PEG donor and a sialyltransferase such as ST3GalIII. In another embodiment the GlcNAc-Asn structures on the peptide/protein backbone can be galactosylated and sialylated as described above, and then further sialylated using CMP-SA-PEG and an  $\alpha$ 2,8-sialyltransferase such as the enzyme encoded by the *Campylobacter jejuni* cst-II gene.

### Herceptin™

#### 32. Glycoconjugation of mithramycin to Herceptin™

This example sets forth the procedures to glycoconjugate a small molecule, such as mithramycin to Fc region glycans of an antibody molecule produced in mammalian cells. Here, the antibody Herceptin™ is used, but one of skill in the art will appreciate that the method can be used with many other antibodies.

**Preparation of Herceptin™-Gal-linker-mithramycin.** Herceptin™ is dissolved at 2.5 mg/mL in 50 mM Tris-HCl, 0.15 M NaCl, 5 mM MnCl<sub>2</sub>, 0.05% NaN<sub>3</sub>, pH 7.2. The solution is incubated with 1 mM UDP-galactose-linker-mithramycin and 0.1 U/mL of galactosyltransferase at 32°C for 2 days to introduce the mithramycin in the Fc region glycans. To monitor the incorporation of galactose, a small aliquot of the reaction has <sup>14</sup>C-galactose-UDP ligand added; the label incorporated into the peptide is separated from the free label by gel filtration on a Toso Haas G3000SW analytical column using PBS buffer (pH 7.1). The radioactive label incorporation into the peptide is quantitated using an in-line radiation detector.

When the reaction is complete, the reaction mixture is purified using a Toso Haas TSK-Gel-3000 preparative column using PBS buffer (pH 7.1) and collecting fractions based on UV absorption. The fractions containing product are combined, concentrated, buffer exchanged and then freeze-dried. The product of the reaction is analyzed using SDS-PAGE and IEF analysis according to the procedures and reagents supplied by Invitrogen. Samples are dialyzed against water and analyzed by MALDI-TOF MS.

Interferon  $\alpha$  and Interferon  $\beta$ 33. GlycoPEGylation of Proteins expressed in Mammalian or Insect Systems:  
EPO, Interferon  $\alpha$  and Interferon  $\beta$ 

5 This example sets forth the preparation of PEGylated peptides that are expressed in mammalian and insect systems.

**Preparation of acceptor from mammalian expression systems.** The peptides to be glycoPEGylated using CMP-sialic acid PEG need to have glycans terminating in galactose. Most peptides from mammalian expression systems will have terminal sialic acid that first  
10 needs to be removed.

**Sialidase digestion.** The peptide is desialylated using a sialidase. A typical procedure involves incubating a 1 mg/mL solution of the peptide in Tris-buffered saline, pH 7.2, with 5 mM  $\text{CaCl}_2$  added, with 0.2 U/mL immobilized sialidase from *Vibrio cholera* (Calbiochem) at 32°C for 24 hours. Microbial growth can be halted either by sterile filtration  
15 or the inclusion of 0.02% sodium azide. The resin is then removed by centrifugation or filtration, and then washed to recover entrapped peptide. At this point, EDTA may be added to the solution to inhibit any sialidase that has leached from the resin.

**Preparation from insect expression systems.** EPO, interferon-alpha, and interferon-beta may also be expressed in non-mammalian systems such as yeast, plants, or  
20 insect cells. The peptides to be glycoPEGylated using CMP-sialic acid PEG need to have glycans terminating in galactose. The majority of the N-glycans on peptides expressed in insect cells, for example, are the trimannosyl core. These glycans are first built out to glycans terminating in galactose before they are acceptors for sialyltransferase.

**Building acceptor glycans from trimannosyl core.** Peptide (1 mg/mL) in Tris-  
25 buffered saline, pH 7.2, containing 5 mM  $\text{MnCl}_2$ , 5 mM UDP-glcNAc, 0.05 U/mL GLCNACT I, 0.05 U/mL GLCNACT II, is incubated at 32°C for 24 hours or until the reaction is substantially complete. Microbial growth can be halted either by sterile filtration or the inclusion of 0.02% sodium azide. After buffer exchange to remove UDP and other small molecules, UDP-galactose and  $\text{MnCl}_2$  are each added to 5 mM, galactosyltransferase is  
30 added to 0.05 U/mL, and is incubated at 32°C for 24H or until the reaction is substantially



complete. Microbial growth can be halted either by sterile filtration or the inclusion of 0.02% sodium azide. The peptides are then ready for glycoPEGylation.

**Building O-linked glycans.** A similar strategy may be employed for interferon alpha to produce enzymatically the desired O-glycan Gal-GalNAc. If necessary, GalNAc linked to serine or threonine can be added to the peptide using appropriate peptide GalNAc transferases (e.g. GalNAc T1, GalNAc T2, T3, T4, etc. ) and UDP-GalNAc. Also, if needed, galactose can be added using galactosyltransferase and UDP-galactose.

**GlycoPEGylation using sialyltransferase.** The glycopeptides (1 mg/mL) bearing terminal galactose in Tris buffered saline + 0.02% sodium azide are incubated with CMP-SA-PEG (0.75 mM) and 0.4 U/mL sialyltransferase (ST3Gal3 or ST3Gal4 for N-glycans on EPO and interferon beta; ST3Gal4, or ST3Gal1 for O-glycans on interferon alpha) at 32°C for 24 hours. Other transferases that may work include the 2,6 sialyltransferase from *Photobacterium damsella*. The acceptor peptide concentration is most preferably in the range of 0.1 mg/mL up to the solubility limit of the peptide. The concentration of CMP-SA-PEG should be sufficient for there to be excess over the available sites, but not so high as to cause peptide solubility problems due to the PEG, and may range from 50 µM up to 5 mM, and the temperature may range from 2°C up to 40°C. The time required for complete reaction will depend on the temperature, the relative amounts of enzyme to acceptor substrate, the donor substrate concentration, and the pH.

#### 34. GlycoPEGylation of Interferon $\alpha$ produced in CHO cells

**Preparation of Asialo-Interferon  $\alpha$ .** Interferon alpha produced from CHO cells is dissolved at 2.5 mg/mL in 50 mM Tris 50 mM Tris-HCl pH 7.4, 0.15 M NaCl, 5 mM CaCl<sub>2</sub> and concentrated to 500 µL in a Centricon Plus 20 centrifugal filter. The solution is incubated with 300 mU/mL Neuraminidase II (*Vibrio cholerae*) for 16 hours at 32 °C. To monitor the reaction a small aliquot of the reaction is diluted with the appropriate buffer and a IEF gel performed. The reaction mixture is then added to prewashed N-(*p*-aminophenyl)oxamic acid-agarose conjugate (800 µL/mL reaction volume) and the washed beads gently rotated for 24 hours at 4 °C. The mixture is centrifuged at 10,000 rpm and the supernatant was collected. The beads are washed 3 times with Tris-EDTA buffer, once with 0.4 mL Tris-EDTA buffer and once with 0.2 mL of the Tris-EDTA buffer and all

supernatants were pooled. The supernatant is dialyzed at 4 °C against 50 mM Tris -HCl pH 7.4, 1 M NaCl, 0.05% NaN<sub>3</sub> and then twice more against 50 mM Tris -HCl pH 7.4, 1 M NaCl, 0.05% NaN<sub>3</sub>. The dialyzed solution is then concentrated using a Centricon Plus 20 centrifugal filter and stored at -20 °C. The conditions for the IEF gel are run according to the procedures and reagents provided by Invitrogen. Samples of native and desialylated G-CSF are dialyzed against water and analyzed by MALDI-TOF MS.

**Preparation of Interferon-alpha-(alpha2,3)-Sialyl-PEG.** Desialylated interferon-alpha is dissolved at 2.5 mg/mL in 50 mM Tris-HCl, 0.15 M NaCl, 0.05% NaN<sub>3</sub>, pH 7.2. The solution is incubated with 1 mM CMP-sialic acid-PEG and 0.1 U/mL of ST3Gal1 at 32°C for 2 days. To monitor the incorporation of sialic acid-PEG, a small aliquot of the reaction had CMP-SA-PEG-fluorescent ligand added; the label incorporated into the peptide is separated from the free label by gel filtration on a Toso Haas G3000SW analytical column using PBS buffer (pH 7.1). The fluorescent label incorporation into the peptide is quantitated using an in-line fluorescent detector. After 2 days, the reaction mixture is purified using Toso Haas G3000SW preparative column using PBS buffer (pH 7.1) and collecting fractions based on UV absorption. The product of the reaction is analyzed using SDS-PAGE and IEF analysis according to the procedures and reagents supplied by Invitrogen. Samples of native and desialylated Interferon-alpha are dialyzed against water and analyzed by MALDI-TOF MS.

**Preparation of Interferon-alpha-(alpha2,8)-Sialyl-PEG.** Interferon-alpha produced in CHO, which contains an alpha2,3-sialylated O-linked glycan, is dissolved at 2.5 mg/mL in 50 mM Tris-HCl, 0.15 M NaCl, 0.05% NaN<sub>3</sub>, pH 7.2. The solution is incubated with 1 mM CMP-sialic acid-PEG and 0.1 U/mL of CST-II at 32°C for 2 days. To monitor the incorporation of sialic acid-PEG, a small aliquot of the reaction has CMP-SA-PEG-fluorescent ligand added; the label incorporated into the peptide is separated from the free label by gel filtration on a Toso Haas G3000SW analytical column using PBS buffer (pH 7.1). The fluorescent label incorporation into the peptide is quantitated using an in-line fluorescent detector. After 2 days, the reaction mixture is purified using a Toso Haas G3000SW preparative column using PBS buffer (pH 7.1) and collecting fractions based on UV absorption. The product of the reaction is analyzed using SDS-PAGE and IEF analysis

according to the procedures and reagents supplied by Invitrogen. Samples of native and PEGylated interferon-alpha are dialyzed against water and analyzed by MALDI-TOF MS.

**Preparation of Interferon-alpha-(alpha2,6)-Sialyl-PEG.** Interferon-alpha, containing only O-linked GalNAc, was dissolved at 2.5 mg/mL in 50 mM Tris-HCl, 0.15 M NaCl, 0.05% NaN<sub>3</sub>, pH 7.2. The solution is incubated with 1 mM CMP-sialic acid-PEG and 0.1 U/mL of ST6GalNAcI or II at 32°C for 2 days. To monitor the incorporation of sialic acid-PEG, a small aliquot of the reaction had CMP-SA-PEG-fluorescent ligand added; the label incorporated into the peptide is separated from the free label by gel filtration on a Toso Haas G3000SW analytical column using PBS buffer (pH 7.1). The fluorescent label incorporation into the peptide is quantitated using an in-line fluorescent detector. After 2 days, the reaction mixture is purified using a Toso Haas G3000SW preparative column using PBS buffer (pH 7.1) and collecting fractions based on UV absorption. The product of the reaction is analyzed using SDS-PAGE and IEF analysis according to the procedures and reagents supplied by Invitrogen. Samples of native and PEGylated interferon-alpha are dialyzed against water and analyzed by MALDI-TOF MS.

### 35. GlycoPEGylation of Interferon-β-1a with PEG (10 kDa) and PEG (20 kDa)

This example illustrates a procedure PEGylate Interferon-β with either PEG (10 kDa) or PEG (20 kDa).

Briefly, Interferon-β-1a (INF-β) was obtained from Biogen (Avonex™). The IFN-β was first purified by Superdex-75 chromatography. The IFN-β was then desialylated with *Vibrio cholerae* sialidase. The INF-β was then PEGylated with SA-PEG (10 kDa) or SA-PEG (20 kDa) and purified with Superdex-200 chromatography.

**Superdex-75 chromatography purification.** INF-β (150 μg) was applied to a Superdex-75 column (Amersham Biosciences, Arlington Heights, IL) and eluted with PBS with 0.5 M NaCl, 0.02 Tween-20, 20 mM histidine and 10% glycerol. The eluant was monitored for absorbance at 280 nm (Figure 172A and 172B) and fractions were collected. Peaks 4 and 5 were pooled, concentrated in an Amicon Ultra 15 spin filter (Millipore, Billerica, MA), and the buffer was exchanged to TBS with 5 mM CaCl<sub>2</sub>, 0.02% Tween-20, 20 mM histidine and 10% glycerol.

**Sialidase Reaction.** The INF- $\beta$  was then desialylated with *Vibrio cholera* sialidase (70 mU/ml, CALBIOCHEM®, EMD Biosciences, Inc., San Diego, CA) on agarose in TBS with 5 mM CaCl<sub>2</sub>, 0.02% Tween-20, 20 mM histidine and 10% glycerol. The reaction was carried out at 32°C for 18 hours. The INF- $\beta$  was removed from the agarose with a 0.22  $\mu$ m Spin-X™ filter (Corning Technology, Inc., Norcross, GA). Figure 173A depicts the MALDI analysis of glycans released from native INF- $\beta$ . The native INF- $\beta$  has many glycoforms containing terminal sialic acid moieties. Figure 173B depicts the MALDI analysis of glycans released from desialylated INF- $\beta$ . The desialylated INF- $\beta$  has primarily one glycoform which is bi-antennary with terminal galactose moieties.

**Lectin Dot-Blot Analysis of Sialylation.** Samples of the INF- $\beta$  from the desialidase reaction were dot-blotted onto nitrocellulose and then blocked with Tris buffered saline (TBS: 0.05M Tris, 0.15M NaCl, pH 7.5) and DIG kit (glycan differentiation kit available from Roche #1 210 238) blocking buffer. Some of the blots were incubated with *Maackia amurensis* agglutinin (MAA) labeled with digoxigenin (DIG) (Roche Applied Science, Indianapolis, IL) to detect  $\alpha$ 2,3-sialylation of INF- $\beta$ . These blots were washed with TBS then incubated with anti-digoxigenin antibody labeled with alkaline phosphatase, then washed again with TBS and developed with NBT/X-phosphate solution, wherein NBT is 4-nitro blue tetrazolium chloride and X-phosphate is 5-bromo-4-chloro-3-indoyl phosphate. The left side of Figure 174 depicts the results of the MAA blot of INF- $\beta$  after the desialylation reaction. The INF- $\beta$  is partially desialylated, as indicated by the decrease in dot development as compared to native INF- $\beta$  in the desialylated samples.

Other blots were incubated with *Erthrina cristagalli* lectin (ECL) labeled with biotin (Vector Laboratories, Burlingame, CA) to detect exposed galactose residues on INF- $\beta$ . After incubation with 2.5  $\mu$ g/ml ECL, the blots were washed in TBS and incubated with streptavidin labeled with alkaline phosphatase. The blots were then washed again and developed. The right side of Figure 174 depicts the ECL blot after development. The increased intensity of the dot of desialylated INF- $\beta$  as compared to the native INF- $\beta$  indicate more exposed galactose moieties and therefore extensive desialylation.

**PEGylation of Desialylated INF- $\beta$  with SA-PEG (10 kDa).** Desialylated INF- $\beta$  (0.05 mg/ml) was PEGylated with ST3Gal3 (50 mU/ml) and CMP-SA-PEG (10 kDa) (250

μM) in an appropriate buffer of TBS + 5 mM CaCl<sub>2</sub>, 0.02% Tween 20, 20 mM histidine, 10% glycerol for 50 hours at 32°C. Figure 175 depicts the SDS-PAGE analysis of the reaction products showing PEGylated INF-β at approximately 98 kDa.

**PEGylation of Desialylated INF-β with SA-PEG (20 kDa).** Desialylated INF-β (0.5 mg/ml) was PEGylated with ST3Gal3 (170 mU/ml) and CMP-SA-PEG (20 kDa) in an appropriate buffer of TBS + 5 mM CaCl<sub>2</sub>, 0.02% Tween 20, 20 mM histidine, 10% glycerol for 50 hours at 32°C. Figure 176 depicts the SDS-PAGE analysis the products of the PEGylation reaction. The PEGylated INF-β has many higher molecular weight bands not found in the unmodified INF-β indicating extensive PEGylation.

**Superdex-200 Purification of INF-β PEGylated with PEG (10 kDa).** The products of the PEGylation reaction were separated on a Superdex-200 column (Amersham Biosciences, Arlington Heights, IL) in PBS with 0.5 NaCl, 0.02 Tween-20, 20 mM histidine and 10% glycerol at 1ml/min and 30 cm/hr flow. The eluant was monitored for absorbance at 280 nm (Figure 177) and fractions were collected. Peaks 3 and 4 were pooled and concentrated in an Amicon Ultra 15 spin filter.

**Bioassay of INF-β PEGylated with PEG (10 kDa).**

The test is inhibition of the proliferation of the lung carcinoma cell line, A549. The A549 cell line are lung carcinoma adherent cells growing in RPMI + 10% FBS at 37°C 5% CO<sub>2</sub>. They can be obtained from ATCC # CCL-185. Wash the cells with 10 ml of PBS and remove the PBS. Add 5 ml of trypsin, incubate for 5 minutes at room temperature or 2 minutes at 37°C. When the cells are detached resuspend into 25 ml of media and count the cells. Dilute the cells at a concentration of 10000 cells/ml and add 200 ul / well (96 wells plate). Incubate for 4 hours at 37°C 5% CO<sub>2</sub>. Prepare 1 ml of IFN B at a concentration of 0.1 ug/ml. Filter it under the hood with a 0.2 um filter. Add 100 ul per well (8 replicates = 1 lane). Incubate for 3 days (do not let the cells go to confluence). Remove 200 ul of media (only 100ul per well left). Add 25 μl of MTT (Sigma) (5 mg/ml filtered 0.22μm). Incubate for 4 hours at 37°C and 5% CO<sub>2</sub>. Aspirate the media gently and add 100 μl of a mixture of isopropanol (100 ml and 6N HCl. Aspirate up and down to homogenize the crystal violet. Read OD 570nm (remove the background at 630 or 690 nm).

Figure 178 depicts the results of the bioassay of the peaks containing INF- $\beta$  PEGylated with PEG (10 kDa) as eluted from the Superdex-200 column.

**Superdex-200 Purification of INF- $\beta$  PEGylated with PEG (20 kDa).** The products of the PEG (20 kDa) PEGylation reaction were separated on a Superdex-200 column (Amersham Biosciences, Arlington Heights, IL) in PBS with 0.5 NaCl, 0.02 Tween-20, 20 mM histidine and 10% glycerol at 1 ml/min flow. The eluant was monitored for absorbance at 280 nm (Figure 179) and fractions were collected. Peak 3 contained most of the INF- $\beta$  PEGylated with PEG (20 kDa).

**Endotoxin test of INF- $\beta$  PEGylated with PEG (20 kDa).**

Limulus Lysate Test was performed, BioWhittaker # 50-647U

**Table 24.** Results of the endotoxin test of INF- $\beta$  PEGylated with PEG (20 kDa).

|                                | Concentration |            |                   |
|--------------------------------|---------------|------------|-------------------|
| INF- $\beta$ with PEG (20 kDa) | 10 EU/ml      | 0.06 mg/ml | 0.16 EU/ $\mu$ g  |
| INF- $\beta$ with PEG (20 kDa) | 1 EU/ml       | 0.07 mg/ml | 0.014 EU/ $\mu$ g |
| Native INF- $\beta$            | 40 EU/ml      | 0.1 mg/ml  | 0.4 EU/ $\mu$ g   |

Remicade<sup>TM</sup>

36. GlycoPEGylation of Remicade<sup>TM</sup> antibody

This example sets forth the procedure to glycoPEGylate a recombinant antibody molecule by introducing PEG molecules to the Fc region glycans. Here Remicade<sup>TM</sup>, a TNF-R:IgG Fc region fusion protein, is the exemplary peptide.

**Preparation of Remicade<sup>TM</sup>-Gal-PEG (10 kDa).** Remicade<sup>TM</sup> is dissolved at 2.5 mg/mL in 50 mM Tris-HCl, 0.15 M NaCl, 5 mM MnCl<sub>2</sub>, 0.05% NaN<sub>3</sub>, pH 7.2. The solution is incubated with 1 mM UDP-galactose-PEG (10 kDa) and 0.1 U/mL of galactosyltransferase at 32°C for 2 days to introduce the PEG in the Fc region glycans. To monitor the incorporation of galactose, a small aliquot of the reaction has <sup>14</sup>C-galactose-UDP ligand added; the label incorporated into the peptide is separated from the free label by gel filtration on a Toso Haas G3000SW analytical column using PBS buffer (pH 7.1). The radioactive label incorporation into the peptide is quantitated using an in-line radiation detector.

When the reaction is complete, the reaction mixture is purified using a Toso Haas TSK-Gel-3000 preparative column using PBS buffer (pH 7.1) and collecting fractions based on UV absorption. The fractions containing product are combined, concentrated, buffer exchanged and then freeze-dried. The product of the reaction is analyzed using SDS-PAGE and IEF analysis according to the procedures and reagents supplied by Invitrogen. Samples are dialyzed against water and analyzed by MALDI-TOF MS.

### Rituxan™

#### 37. Glycoconjugation of geldanamycin to Rituxan™

This example sets forth the glycoconjugation of a small molecule, such as geldanamycin, to the Fc region glycans of an antibody produced in CHO cells, such as Rituxan™. Here, the antibody Rituxan™ is used, but one of skill in the art will appreciate that the method can be used with many other antibodies.

**Preparation of Rituxan™-Gal-linker-geldanamycin.** Rituxan™ is dissolved at 2.5 mg/mL in 50 mM Tris-HCl, 0.15 M NaCl, 5 mM MnCl<sub>2</sub>, 0.05% NaN<sub>3</sub>, pH 7.2. The solution is incubated with 1 mM UDP-galactose-linker-geldanamycin and 0.1 U/mL of galactosyltransferase at 32°C for 2 days to introduce the geldanamycin in the Fc region glycans. To monitor the incorporation of galactose, a small aliquot of the reaction has <sup>14</sup>C-galactose-UDP ligand added; the label incorporated into the peptide is separated from the free label by gel filtration on a Toso Haas G3000SW analytical column using PBS buffer (pH 7.1). The radioactive label incorporation into the peptide is quantitated using an in-line radiation detector.

When the reaction is complete, the reaction mixture is purified using a Toso Haas TSK-Gel-3000 preparative column using PBS buffer (pH 7.1) and collecting fractions based on UV absorption. The fractions containing product are combined, concentrated, buffer exchanged and then freeze-dried. The product of the reaction is analyzed using SDS-PAGE and IEF analysis according to the procedures and reagents supplied by Invitrogen. Samples are dialyzed against water and analyzed by MALDI-TOF MS.

Rnase38. Remodeling high mannose N-glycans to hybrid and complex N-glycans:  
Bovine pancreatic RNase

This example sets forth the preparation of bovine pancreas RNase with hybrid or  
5 complex N-glycans. The high mannose N-linked glycans of the RNase are enzymatically  
digested and elaborated to create hybrid N-linked glycans. Additionally, the high mannose  
N-linked glycans of the RNase are enzymatically digested and elaborated to create complex  
N-linked glycans.

High mannose structures of *N*-linked oligosaccharides in glycopeptides can be  
10 modified to hybrid or complex forms using the combination of  $\alpha$ -mannosidases and  
glycosyltransferases. This example summarizes the results in such efforts using a simple *N*-  
Glycan as a model substrate.

Ribonuclease B (RNaseB) purified from bovine pancreas (Sigma) is a glycopeptide  
consisting of 124 amino acid residues. It has a single potential *N*-glycosylation site modified  
15 with high mannose structures. Due to its simplicity and low molecular weight (13.7 kDa to  
15.5 kDa), ribonuclease B is a good candidate to demonstrate the feasibility of the *N*-Glycan  
remodeling from high mannose structures to hybrid or complex *N*-linked oligosaccharides.  
The MALDI-TOF spectrum of RNaseB (Figure 180A) and HPLC profile for the  
oligosaccharides cleaved from RNaseB by *N*-Glycanase (Figure 180B) indicated that, other  
20 than a small portion of the non-modified peptide, the majority of *N*-glycosylation sites of the  
peptide are modified with high mannose oligosaccharides consisting of 5 to 9 mannose  
residues.

**Conversion of high mannose N-Glycans to hybrid N-Glycans.** High mannose *N*-  
Glycans were converted to hybrid *N*-Glycans using the combination of  $\alpha$ 1,2-mannosidase,  
25 GlcNAcT-I ( $\beta$ -1,2-*N*-acetyl glucosaminyl transferase), GalT-I ( $\beta$ 1,4-galactosyltransferase) and  
 $\alpha$ 2,3-sialyltransferase /or  $\alpha$ 2,6-sialyltransferase as shown in Figure 181.

As an example, high mannose structures in RNaseB were successfully converted to  
hybrid structures.

Man<sub>5</sub>GlcNAc<sub>2</sub>-R was obtained from Man<sub>5-9</sub>GlcNAc<sub>2</sub>-R catalyzed by a single  $\alpha$ 1,2-  
30 mannosidase cloned from *Trichoderma reesei* (Figure 182). RNase B (1 g, about 67  $\mu$ mol)  
was incubated at 30°C for 45 hr with 15 mU of the recombinant *T. reesei*  $\alpha$ 1,2-mannosidase



in MES buffer (50 mM, pH 6.5) in a total volume of 10 mL. Man<sub>6-9</sub>GlcNAc<sub>2</sub>-protein structures have been successfully converted to Man<sub>5</sub>GlcNAc<sub>2</sub>-protein with high efficiency by the recombinant mannosidase.

Alternately, Man<sub>5</sub>GlcNAc<sub>2</sub>-R was obtained from Man<sub>5-9</sub>GlcNAc<sub>2</sub>-R catalyzed by a single  $\alpha$ 1,2-mannosidase purified from *Aspergillus saitoi* (Figure 183). RNase B (40  $\mu$ g, about 2.7 nmol) was incubated at 37°C for 42.5 hr with 25  $\mu$ U of the commercial *A. saitoi*  $\alpha$ 1,2-mannosidase (Glyko or CalBioChem) in NaOAc buffer (100 mM, pH 5.0) in a total volume of 20  $\mu$ l. Man<sub>6-9</sub>GlcNAc<sub>2</sub>-protein structures were successfully converted to Man<sub>5</sub>GlcNAc<sub>2</sub>-protein by the commercially available mannosidase. However, a new peak corresponding to the GlcNAc-protein appears in the spectrum, indicating the possible contamination of endoglycosidase H in the preparation. Although several mammalian alpha-mannosidases were required to achieve this step, the fungal  $\alpha$ 1,2-mannosidase was very efficient to remove all  $\alpha$ 1,2-linked mannose residues.

GlcNAcT-I then added a GlcNAc residue to the Man<sub>5</sub>GlcNAc<sub>2</sub>-R (Figure 184). The reaction mixture after the *T. reesei*  $\alpha$ 1,2-mannosidase reaction containing RNase B (600  $\mu$ g, about 40 nmol) was incubated with non-purified recombinant GlcNAcT-I (34 mU) in MES buffer (50 mM, pH 6.5) containing MnCl<sub>2</sub> (20 mM) and UDP-GlcNAc (5 mM) in a total volume of 400  $\mu$ l. at 37°C for 42 hr. A GlcNAc residue was quantitatively added to Man<sub>5</sub>GlcNAc<sub>2</sub>-protein by the recombinant GlcNAcT-I.

A Gal residue was then added using GalT 1 (Figure 185). The reaction mixture after the GnT-I reaction containing RNase B (120  $\mu$ g, about 8 nmol) was incubated at 37°C for 20 hr with 3.3 mU of the recombinant GalT-1 in Tris-HCl buffer (100 mM, pH 7.3) containing UDP-Gal (7.5 mM) and MnCl<sub>2</sub> (20 mM) in a total volume of 100  $\mu$ l. A Gal residue was added to about 98% of the GlcNAc-Man<sub>5</sub>GlcNAc<sub>2</sub>-protein by the recombinant GalT 1.

The next step was the addition of a sialic acid using an  $\alpha$ 2,3-sialyltransferase or an  $\alpha$ 2,6-sialyltransferase (Figure 186). As an example, ST3Gal III, an  $\alpha$ 2,3-sialyltransferase was used. The reaction mixture after the GalT-1 reaction containing RNase B (13  $\mu$ g, about 0.87 nmol) was incubated at 37°C for 16 hr with 8.9 mU of recombinant ST3Gal III in Tris-HCl buffer (100 mM, pH 7.3) containing CMP-Sialic acid (5 mM) and MnCl<sub>2</sub> (20 mM) in a total volume of 20  $\mu$ l. A sialic acid residue was added to about 90% of the Gal-GlcNAc-

Man<sub>5</sub>GlcNAc<sub>2</sub>-protein by recombinant ST3Gal III using CMP-SA as the donor. The yield can be further improved by adjusting the reaction conditions.

For convenience, no purification or dialysis step was required after each reaction described above. More interesting, GalT 1 and ST3Gal III can be combined in a one-pot reaction. Similar yields were obtained as compared with the separate reactions. The reaction mixture after the GlcNAcT-I reaction containing RNase B (60 µg, about 4 nmol) was incubated at 37°C for 20 hr with 1.7 mU of recombinant GalT 1, 9.8 mU of recombinant ST3Gal III in Tris-HCl buffer (100 mM, pH 7.3) containing UDP-Gal (7.5 mM), CMP-sialic acid (5 mM) and MnCl<sub>2</sub> (20 mM) in a total volume of 60 µl.

As shown in Figure 187, SA-PEG (10 kDa) was successfully added to the RNaseB. The reaction mixture after the GalT-1 reaction containing RNase B (6.7 µg, about 0.45 nmol) was dialyzed against H<sub>2</sub>O for 1 hour at room temperature and incubated at 37°C for 15.5 hours with 55 mU of the recombinant ST3Gal III in Tris-HCl buffer (50 mM, pH 7.3) containing CMP-SA-PEG (10 kDa) (0.25 mM) and MnCl<sub>2</sub> (20 mM) in a total volume of 20 µl. PEG-modified sialic acid residues were successfully added to the Gal-GlcNAc-Man<sub>5</sub>GlcNAc<sub>2</sub>-peptide by the recombinant ST3Gal III. The yield can be further improved by adjusting the reaction conditions.

**Conversion of high mannose N-Glycans to complex N-Glycans.** To achieve this conversion, a GlcNAcβ1,2Man<sub>3</sub>GlcNAc<sub>2</sub>-peptide intermediate is obtained. As shown in Figure 188, there are at least four feasible routes to carry out the reaction from Man<sub>5</sub>GlcNAc<sub>2</sub>-peptide to this intermediate:

**Route I:** The Man<sub>5</sub>GlcNAc<sub>2</sub>-peptide produced by the fungal α1,2 mannosidase is a substrate of GlcNAc transferase I (GlcNAcT-I, enzyme 2) which adds one GlcNAc. The terminal α1,3- and α1,6-linked mannose residues of GlcNAcMan<sub>5</sub>GlcNAc<sub>2</sub>-peptide is removed by Golgi α-mannosidase II (ManII, enzyme 5). This route is a part of the natural pathway for the processing of N-linked oligosaccharides carried out in higher organisms.

**Route II:** Two mannose residues are first removed by an α-mannosidase (enzyme 6), then a GlcNAc is added by GlcNAcT-I (enzyme 2). Other than its natural acceptor Man<sub>5</sub>GlcNAc<sub>2</sub>-R, GlcNAcT-I can also recognize Man<sub>3</sub>GlcNAc<sub>2</sub>-R as its substrate and add one GlcNAc to the mannose core structure to form GlcNAcMan<sub>3</sub>GlcNAc<sub>2</sub>-peptide.

**Route III:** The  $\alpha$ 1,6-linked mannose is removed by an  $\alpha$ 1,6-mannosidase, followed by the addition of GlcNAc by GlcNAcT-I and removal of the terminal  $\alpha$ 1,3-linked mannose by an  $\alpha$ 1,3-mannosidase. From the experimental data obtained, GlcNAcT-I can recognize this Man<sub>4</sub>GlcNAc<sub>2</sub>-peptide as acceptor and add one GlcNAc residue to form

5 GlcNAcMan<sub>4</sub>GlcNAc<sub>2</sub>-peptide.

**Route IV:** Similar to Route III,  $\alpha$ 1,3-linked mannose is removed by an  $\alpha$ 1,3-mannosidase, followed by GlcNAcT-I reaction. Then the terminal  $\alpha$ 1,6-linked mannose can be removed by an  $\alpha$ 1,6-mannosidase.

After the function of GlcNAcT-I (responsible for the addition of the GlcNAc  $\beta$ 1,2-linked to the  $\alpha$ 1,3-mannose on the mannose core) and GlcNAcT-II (responsible for the addition of a second GlcNAc  $\beta$ 1,2-linked to the  $\alpha$ 1,6-mannose on the mannose core), the GlcNAc<sub>2</sub>Man<sub>3</sub>GlcNAc<sub>2</sub>-peptide can be processed by GalT 1 and sialyltransferase to form bi-antennary complex N- Glycans. Other GlcNAc transferases such as GlcNAcT-IV, GlcNAcT-V, and/or GlcNAcT-VI (Figure 188 and Figure 189) can also glycosylate the

10 linked to the  $\alpha$ 1,3-mannose on the mannose core) and GlcNAcT-II (responsible for the addition of a second GlcNAc  $\beta$ 1,2-linked to the  $\alpha$ 1,6-mannose on the mannose core), the GlcNAc<sub>2</sub>Man<sub>3</sub>GlcNAc<sub>2</sub>-peptide can be processed by GalT 1 and sialyltransferase to form bi-antennary complex N- Glycans. Other GlcNAc transferases such as GlcNAcT-IV, GlcNAcT-V, and/or GlcNAcT-VI (Figure 188 and Figure 189) can also glycosylate the

15 GlcNAc<sub>2</sub>Man<sub>3</sub>GlcNAc<sub>2</sub>-peptide. Additional glycosylation by the GalT 1 and sialyltransferases will form multi-antennary complex N-glycans. The enzyme GlcNAcT-III catalyzes the insertion of a bisecting GlcNAc, thus preventing the actions of ManII and subsequent action of transferases GlcNAcT-II, GlcNAcT-IV and GlcNAcT-V.

## 20 Tissue-Type Plasminogen Activator (TPA)

### 39. Fucosylation of TPA to create Sialyl Lewis X

This example sets forth the preparation of Tissue Tissue-type Plasminogen Activator (TPA) with N-linked sialyl Lewis X antigen.

**Sialylation.** TPA expressed in mammalian cells will often contain a majority of the glycans terminating in sialic acid, but to ensure complete sialylation, it would be beneficial to first perform an *in vitro* sialylation. TPA in a suitable buffer (most preferably between pH 5.5 and 9, for example Tris buffered saline, pH 7.2) is incubated with CMP sialic acid and sialyltransferase for a time sufficient to convert any glycans lacking sialic acid to sialylated species. Typical conditions would be 1 mg/mL TPA, 3 mM CMP sialic acid, 0.02 U/mL

25 glycans terminating in sialic acid, but to ensure complete sialylation, it would be beneficial to first perform an *in vitro* sialylation. TPA in a suitable buffer (most preferably between pH 5.5 and 9, for example Tris buffered saline, pH 7.2) is incubated with CMP sialic acid and sialyltransferase for a time sufficient to convert any glycans lacking sialic acid to sialylated species. Typical conditions would be 1 mg/mL TPA, 3 mM CMP sialic acid, 0.02 U/mL

30 ST3Gal3, 32°C for 24 hours. Microbial growth can be halted either by sterile filtration or the

inclusion of 0.02% sodium azide. The TPA concentration is most preferably in the range 0.1 mg/mL up to the solubility limit of the peptide. The concentration of CMP-SA should be sufficient for there to be excess over the available sites, and might range from 50  $\mu$ M up to 50 mM, and the temperature from 2°C up to 40°C. The time required for complete reaction will depend on the temperature, the relative amounts of enzyme to acceptor substrate, the donor substrate concentration, and the pH. Other sialyltransferases that may be capable of adding sialic acid in 2,3 linkage include ST3Gal4; microbial transferases could also be used.

**Fucosylation.** Typical conditions for fucosylation would be 1 mg/mL TPA, 3 mM GDP-fucose, 0.02 U/mL FTVI, 5 mM MnCl<sub>2</sub>, 32°C for 24H in Tris buffered saline.

Microbial growth can be halted either by sterile filtration or the inclusion of 0.02% sodium azide. The TPA concentration is most preferably in the range 0.1 mg/mL up to the solubility limit of the peptide. The concentration of GDP-fucose should be sufficient for there to be excess over the available sites, and might range from 50  $\mu$ M up to 50 mM, and the temperature from 2°C up to 40°C. The time required for complete reaction will depend on the temperature, the relative amounts of enzyme to acceptor substrate, the donor substrate concentration, and the pH. Other fucosyltransferases that may be capable of making sialyl Lewis x include FTVII, FTV, FTIII, as well as microbial transferases could also be used.

#### 40. Trimming of high mannose to tri-mannose core structure: Tissue-type Plasminogen Activator produced in CHO

This example sets forth the preparation of Tissue-type Plasminogen Activator with a trimannose core by trimming back from a high mannose glycan.

Tissue-type plasminogen activator (TPA) is currently produced in Chinese Hamster Ovary (CHO) cells and contains a low amount of high mannose N-linked oligosaccharide.

The mannoses can be trimmed down using a variety of the specific mannosidases. The first step is to generate Man5GlcNAc2(Fuc0-1) from Man9GlcNAc2(Fuc0-1). This can be done using mannosidase I. Then either GlcNAcT1 (GlcNAc transferase I) is used to make GlcNAc1Man5GlcNAc2(Fuc0-1) or Mannosidase III is used to make Man3GlcNAc2(Fuc0-1). From Man3GlcNAc2(Fuc0-1), GlcNAc1Man3GlcNAc2(Fuc0-1) can be produced using GlcNAcT1 or from GlcNAc1Man5GlcNAc2(Fuc0-1), GlcNAc1Man3GlcNAc2(Fuc0-1) can be produced using Mannosidase II. GlcNAc1Man3GlcNAc2(Fuc0-1) is then converted into

GlcNAc2Man3GlcNAc2(Fuc0-1) using GlcNAcTransferase II (GlcNAcTII). The two terminal GlcNAc residues are then galactosylated using GalTI and then sialylated with SA-PEG using ST3GalIII.

Conversely, TPA can be produce in yeast or fungal systems. Similar processing  
5 would be required for fungal derived material.

#### 41. Generation and PEGylation of GlcNAc-Asn structures: TPA produced in Yeast

This example sets forth the preparation of PEGylated GlcNAc-Asn structures on a  
10 peptide such as TPA expressed in yeast.

Yeast expression is expected to result in a TPA which contains a single N-linked mannan-type structure. This recombinant glycoprotein is first treated with endoglycosidase H to generate GlcNAc structures on the asparagine (Asn) residues on the peptide.

The GlcNAc-Asn structures on the peptide/protein backbone are then be modified  
15 with galactose or galactose-PEG using UDP-galactose or UDP-galactose-6-PEG, respectively, and a galactosyltransferase such as GalT1. In one case, the galactose-PEG is the terminal residue. In the second case, the galactose is further modified with SA-PEG using a CMP-SA-PEG donor and a sialyltransferase such as ST3GalIII. In another embodiment, the GlcNAc-Asn structures on the peptide/protein backbone may be  
20 galactosylated and sialylated as described above, and then further sialylated using CMP-SA-PEG and an  $\alpha$ 2,8-sialyltransferase such as the enzyme encoded by the *Campylobacter jejuni* cst-II gene.

#### Transferrin

#### 25 42. GlycoPEGylation of Transferrin

This example sets forth the preparation of asialotransferrin and its sialylation with PEG-CMP-sialic acid.

**Preparation of Asialo-transferrin.** Human-derived holo-Transferrin, (10 mg) was dissolved in 500  $\mu$ L of 50 mM NaOAc, 5 mM CaCl<sub>2</sub>, pH 5.5. To this solution was added  
30 500 mU Neuraminidase II (*Vibrio cholerae*) and the reaction mixture was shaken gently for 20.5 hours at 37 °C. The reaction mixture was added to the prewashed N-(p-

aminophenyl)oxamic acid-agarose conjugate (600  $\mu$ L) and the washed beads gently rotated for 24 hours at 4 °C. The mixture was centrifuged at 10,000 rpm and the supernatant was collected. The reaction mixture was adjusted to 5 mM EDTA by addition of 100  $\mu$ L of 30 mM EDTA to the washed beads, which were gently rotated for 20 hours at 4 °C. The suspension was centrifuged for 2 minutes at 10,000 rpm and the supernatant was collected. The beads were washed 5 times with 0.35 mL of 50 mM NaOAc, 5 mM CaCl<sub>2</sub>, 5 mM EDTA, pH 5.5 and all supernatants were pooled. The enzyme solution was dialyzed twice at 4 °C into 15 mM Tris-HCl, 1 M NaCl, pH 7.4. 0.3 mL of the transferrin solution (3.3 mL total) was removed and dialyzed twice against water. The remainder was dialyzed twice more at 4 °C against phosphate buffered saline. The dialyzed solution was stored at -20 °C. Protein samples were analyzed by IEF Electrophoresis. Samples (9  $\mu$ L, 25  $\mu$ g) were diluted with 16  $\mu$ L Tris buffer and mixed with 25  $\mu$ L of the sample loading buffer and applied to Isoelectric Focusing Gels (pH 3-7). Gels were run and fixed using standard procedures. Gels were stained with Colloidal Blue Stain.

**Sialyl-PEGylation of asialo-Transferrin.** Desialylated transferrin (250  $\mu$ g) and CMP-sialic acid or CMP-SA-PEG (1 kDa or 10 kDa)(0.05  $\mu$ mol) were dissolved in 69  $\mu$ L 50 mM Tris-HCl, 0.15 M NaCl, 0.05% NaN<sub>3</sub>, pH 7.2 in 1.5 mL plastic tubes. The tubes were vortexed briefly and 100 mU ST3Gal3 (90  $\mu$ L) were added (total volume 250  $\mu$ L). The tubes were vortexed again and mixed gently for 24 hours at 32 °C. The reactions were stopped by freezing at -80 °C. Novex Tris-Glycine 8-16% 1 mm gels were used for SDS PAGE analysis (Figure 190). Samples (25  $\mu$ L, 25  $\mu$ g) were mixed with 25  $\mu$ L of sample loading buffer and 0.4  $\mu$ L of  $\beta$ -mercaptoethanol and heated for 6 minutes at 85 °C. Gels were run using standard conditions and stained with Colloidal Blue Stain. IEF gels were also performed as described above Figure 191). Samples were also dialyzed against water analyzed by MALDI-TOF.

**Results.** MALDI was also performed. Native transferrin (78729); asialotransferrin (78197); resialylated transferrin (79626/80703); with SA-PEG 1k (79037 (1); 80961 (2); 82535 (3); 84778 (4)); with SA-PEG 5k (90003 (2); 96117 (3); 96117 (4)); with SA-PEG 10k (100336 (2); 111421 (3); 122510 (4)).

#### 43. Transferrin-GDNF

This example sets forth the procedures for the glycoconjugation of proteins, and in particular, transferrin is glycoconjugated to GDNF. Transferrin-SA-Linker-Gal-UDP is prepared from transferrin. The galactose residue is removed from GDNF glycans, and  
5 Transferrin-SA-Linker-Gal-UDP is conjugated to GDNF glycans using a galactosyltransferase.

**Preparation of agalacto-GDNF.** GDNF produced in NSO cells (NSO murine myeloma cells) is dissolved at 2.5 mg/mL in 50 mM Tris 50 mM Tris-HCl pH 7.4, 0.15 M NaCl, and is incubated with 300 mU/mL beta-galactosidase-agarose conjugate for 16 hours at  
10 32°C. To monitor the reaction a small aliquot of the reaction is diluted with the appropriate buffer and a IEF gel performed according to Invitrogen procedures. The mixture is centrifuged at 10,000 rpm and the supernatant is collected. The supernatant is dialyzed at 4 °C against 50 mM Tris -HCl pH 7.4, 1 M NaCl, 0.05% NaN<sub>3</sub> and then twice more against 50 mM Tris -HCl pH 7.4, 1 M NaCl, 0.05% NaN<sub>3</sub>. The dialyzed solution is then concentrated  
15 using a Centricon Plus 20 centrifugal filter and stored at -20 °C. The conditions for the IEF gel are run according to the procedures and reagents provided by Invitrogen. Samples are dialyzed against water and analyzed by MALDI-TOF MS.

**Preparation of Transferrin-SA-Linker-Gal-UDP.** Asialo-transferrin is dissolved at 2.5 mg/mL in 50 mM Tris-HCl, 0.15 M NaCl, 0.05% NaN<sub>3</sub>, pH 7.2. The solution is  
20 incubated with CMP-sialic acid-linker-Gal-UDP (molar amount to add 1 molar equivalent of nucleotide sugar to transferrin) and 0.1 U/mL of ST3Gal3 at 32°C for 2 days. To monitor the incorporation of sialic acid, a small aliquot of the reaction has <sup>14</sup>C-SA-UDP ligand added; the label incorporated into the peptide is separated from the free label by gel filtration on a Toso Haas G3000SW analytical column using PBS buffer (pH 7.1). The radioactive label  
25 incorporation into the peptide is quantitated using an in-line radiation detector.

The solution is incubated with 5 mM CMP-sialic acid and 0.1 U/mL of ST3Gal3 (to cap any unreacted transferrin glycans) at 32°C for 2 days. The incorporation into the peptide is quantitated using an in-line UV detector. After 2 days, the reaction mixture is purified using a Toso Haas G3000SW preparative column using PBS buffer (pH 7.1) and collecting  
30 fractions based on UV absorption. The product of the reaction is analyzed using SDS-PAGE

and IEF analysis according to the procedures and reagents supplied by Invitrogen. Samples are dialyzed against water and analyzed by MALDI-TOF MS.

**Preparation of Transferrin-SA-Linker-Gal-GDNF.** The transferrin-SA-Linker-Gal-UDP prepared as described above is dissolved at 2.5 mg/mL in 50 mM Tris-HCl, 0.15 M NaCl, 5 mM MnCl<sub>2</sub>, 0.05% NaN<sub>3</sub>, pH 7.2. The solution is incubated with 2.5 mg/mL agalacto-GDNF and 0.1 U/mL of galactosyltransferase at 32°C for 2 days. To monitor the incorporation of galactose, a small aliquot of the reaction has <sup>14</sup>C-galactose-UDP ligand added; the label incorporated into the peptide is separated from the free label by gel filtration on a Toso Haas G3000SW analytical column using PBS buffer (pH 7.1). The radioactive label incorporation into the peptide is quantitated using an in-line radiation detector.

When the reaction is complete, the solution is incubated with 5 mM UDP-Gal and 0.1 U/mL of galactosyltransferase (to cap any unreacted transferrin glycans) at 32°C for 2 days followed by addition of 5 mM CMP-SA and 0.1 U/mL of ST3Gal3. After 2 additional days, the reaction mixture is purified using a Toso Haas G3000SW preparative column using PBS buffer (pH 7.1) collecting fractions based on UV absorption. The product of the reaction is analyzed using SDS-PAGE and IEF analysis according to the procedures and reagents supplied by Invitrogen. Samples are dialyzed against water and analyzed by MALDI-TOF MS.

The disclosures of each and every patent, patent application, and publication cited herein are hereby incorporated herein by reference in their entirety.

While this invention has been disclosed with reference to specific embodiments, it is apparent that other embodiments and variations of this invention may be devised by others skilled in the art without departing from the true spirit and scope of the invention. The appended claims are intended to be construed to include all such embodiments and equivalent variations.



What is claimed:

1. An EPO peptide comprising one or more glycans, having a glycoconjugate molecule covalently attached to said peptide.

2. The EPO peptide of claim 1, wherein said one or more glycans is a  
5 monoantennary glycan.

3. The EPO peptide of claim 1, wherein said one or more glycans is a biantennary glycan.

4. The EPO peptide of claim 1, wherein said one or more glycans is a triantennary glycan.

10 5. The EPO peptide of claim 1, wherein said one or more glycans is at least a triantennary glycan.

6. The EPO peptide of claim 1, wherein said one or more glycans comprises at least two glycans comprising a mixture of mono or multiantennary glycans.

15 7. The EPO peptide of claim 1, wherein said one or more glycans is selected from an N-linked glycan and an O-linked glycan.

8. The EPO peptide of claim 1, wherein said one or more glycans is at least two glycans selected from an N-linked and an O-linked glycan.

9. The EPO peptide of claim 1, wherein said peptide is expressed in a cell selected from the group consisting of a prokaryotic cell and a eukaryotic cell.

20 10. The EPO peptide of claim 9, wherein said eukaryotic cell is selected from the group consisting of a mammalian cell, an insect cell and a fungal cell.

11. The EPO peptide of claim 10, wherein said fungal cell is a yeast cell.

12. A glycoPEGylated EPO peptide comprising an EPO peptide and at least one glycan and at least one poly(ethylene glycol) molecule covalently attached to said glycan,

wherein said poly(ethylene glycol) molecule is added to said EPO peptide using a glycosyltransferase.

13. The glycoPEGylated EPO peptide of claim 12, comprising at least one  
5 mono-antennary glycan.

14. The glycoPEGylated EPO peptide of claim 12, wherein all of said glycans are N-linked and are mono-antennary.

10 15. The glycoPEGylated EPO peptide of claim 12, wherein all of said glycans are N-linked and at least one of said glycans comprise said poly(ethylene glycol).

16. The glycoPEGylated EPO peptide of claim 15, wherein more than one of said glycans comprises said poly(ethylene glycol).  
15

17. The glycoPEGylated EPO peptide of claim 12, wherein all of said glycans are N-linked and all of said glycans comprise said poly(ethylene glycol).

18. The glycoPEGylated EPO peptide of claim 12, comprising at least three  
20 mono-antennary glycans having said poly(ethylene glycol) covalently attached thereto.

19. A glycoPEGylated EPO peptide, wherein said EPO peptide comprises three or more glycans.

25 20. The glycoPEGylated EPO peptide of claim 9, wherein at least one of said glycans comprises said poly(ethylene glycol) covalently attached thereto.

21. The glycoPEGylated EPO peptide of claim 18, wherein more than one of said glycans comprises said poly(ethylene glycol) covalently attached thereto.  
30

22. The glycoPEGylated EPO peptide of claim 18, wherein all of said glycans comprise said poly(ethylene glycol) covalently attached thereto.

5 23. The glycoPEGylated EPO peptide of claim 12 wherein said poly(ethylene glycol) is linked to at least one sugar moiety selected from the group consisting of fucose (Fuc), N-acetylglucosamine (GlcNAc), galactose (Gal) and a sialic acid (SA).

24. The glycoPEGylated EPO peptide of claim 23, wherein said sialic acid is N-acetylneuraminic acid.

10 25. The glycoPEGylated EPO peptide of claim 12, wherein said EPO peptide does not comprise an O-linked glycan.

26. The glycoPEGylated EPO peptide of claim 12 wherein said EPO peptide comprises at least one O-linked glycan.

15

27. The glycoPEGylated EPO peptide of claim 26, wherein said O-linked peptide comprises said poly(ethylene glycol) covalently attached thereto.

20 28. The glycoPEGylated EPO peptide of claim 27, wherein said EPO peptide is recombinantly expressed in a cell.

29. The glycoPEGylated EPO peptide of claim 28, wherein said cell is selected from the group consisting of an insect cell, a fungal cell and a mammalian cell.

25 30. The glycoPEGylated EPO peptide of claim 29, wherein said fungal cell is a yeast cell.

31. The glycoPEGylated EPO peptide of claim 29, wherein said cell is an insect cell.

30

32. The glycoPEGylated EPO peptide of claim 29, wherein said cell is a yeast cell.

33. The glycoPEGylated EPO peptide of claim 29, wherein said cell is a mammalian cell.

34. The glycoPEGylated EPO peptide of claim 33, wherein said mammalian cell is a CHO cell.

35. The glycoPEGylated EPO peptide of claim 12, wherein said poly(ethylene glycol) has a molecular weight selected from the group consisting of about 1 kDa, 2 kDa, 5 kDa, 10 kDa, 20 kDa, 30 kDa and 40 kDa.

36. The glycoPEGylated EPO peptide of claim 35, wherein said poly(ethylene glycol) has a molecular weight of 20 kDa.

37. The glycoPEGylated EPO peptide of claim 12, wherein said EPO peptide is selected from the group consisting of a naturally occurring EPO peptide and a mutated EPO peptide.

38. The glycoPEGylated EPO peptide of claim 37, wherein said mutated EPO peptide comprises the amino acid sequence of SEQ ID NO:73 having at least one mutation selected from the group consisting of Arg<sup>139</sup> to Ala<sup>139</sup>, Arg<sup>143</sup> to Ala<sup>143</sup> and Lys<sup>154</sup> to Ala<sup>154</sup>.

39. A method of making a glycoPEGylated EPO peptide, said method comprising the step of:

(a) contacting an EPO peptide with a mixture comprising a nucleotide sugar covalently linked to poly(ethylene glycol) and a glycosyltransferase under conditions sufficient to transfer said poly(ethylene glycol) to said EPO peptide.

40. The method of claim 39, wherein the sugar of said nucleotide sugar is selected from the group consisting of fucose (Fuc), N-acetylglucosamine (GlcNAc), galactose (Gal) and a sialic acid (SA).

5           41. The method of claim 40, wherein said sialic acid is N-acetylneuraminic acid (NAN).

          42. The method of claim 39, wherein said poly(ethylene glycol) has a molecular weight selected from the group consisting of about 1 kDa, 2 kDa, 5 kDa, 10 kDa,  
10   20 kDa, 30 kDa and 40 kDa.

          43. The method of claim 42, wherein said poly(ethylene glycol) has a molecular weight of 20 kDa.

15           44. The method of claim 39, wherein said EPO peptide is recombinantly expressed in a cell.

          45. The method of claim 44, wherein said cell is selected from the group consisting of an insect cell, a fungal cell and a mammalian cell.

20           46. The method of claim 45, wherein said cell is an insect cell.

          47. The method of claim 45, wherein said cell is a yeast cell.

25           48. The method of claim 45, wherein said cell is a mammalian cell.

          49. The method of claim 48, wherein said mammalian cell is a CHO cell.

30           50. The method of claim 39, wherein said EPO peptide is selected from the group consisting of a naturally occurring EPO peptide and a mutated EPO peptide.

51. The method of claim 50, wherein said mature EPO peptide has the sequence of SEQ ID NO:73.

52. The method of claim 50, wherein said mutated EPO peptide comprises the amino acid sequence of SEQ ID NO: 73 having at least one mutation selected from the group consisting of Arg<sup>139</sup> to Ala<sup>139</sup>, Arg<sup>143</sup> to Ala<sup>143</sup> and Lys<sup>154</sup> to Ala<sup>154</sup>.

53. The method of claim 39, wherein before step (a):

(b) contacting said EPO peptide with a mixture comprising a nucleotide-N-acetylglucosamine (GlcNAc) molecule and an N-acetylglucosamine transferase (GnT) for which the nucleotide-GlcNAc is a substrate under conditions sufficient to form a bond between said GlcNAc and said EPO, wherein said GnT is selected from the group consisting of GnT I, GnT II, GnT III, GnT IV, GnT V and GnT VI.

54. The method of claim 53, wherein said mixture comprises one GnT selected from the group consisting of GnT I, GnT II, GnT IV, GnT V and GnT VI.

55. The method of claim 54, wherein said GnT is GnT I.

56. The method of claim 54, wherein said GnT is GnT II.

57. The method of claim 39, wherein said glycoPEGylated EPO peptide comprises at least one mono-antennary glycan.

58. The method of claim 39, wherein the sugar of said nucleotide sugar is galactose and said glycosyltransferase is galactosyl transferase I (GalT I).

59. The method of claim 53, wherein before step (a) but after step (b):  
(c) contacting said EPO peptide with a mixture comprising a nucleotide galactose (Gal) and galactosyl transferase I (GalT I) under conditions sufficient to transfer galactose to said EPO peptide.

60. The method of claim 39, wherein in step (a), the sugar of said nucleotide sugar is sialic acid and said glycosyltransferase is a sialyltransferase.

5           61. The method of claim 60, wherein said sialic acid is N-acetylneuraminic acid (NAN).

62. The method of claim 60, wherein said sialyltransferase is selected from the group consisting of  $\alpha(2,3)$ sialyltransferase,  $\alpha(2,6)$ sialyltransferase and  
10   (2,8)sialyltransferase.

63. A glycoPEGylated EPO peptide made by the method of claim 39.

64. A glycoPEGylated EPO peptide, said EPO peptide comprising the  
15   sequence of SEQ ID NO:73.

65. A glycoPEGylated EPO peptide, said EPO peptide comprising the sequence of SEQ ID NO:73 and further comprising a mutation in said sequence.

20           66. A method of making a glycoPEGylated EPO peptide, said method comprising the steps of:

          (a) contacting an EPO peptide with a mixture comprising a nucleotide sugar covalently linked to poly(ethylene glycol) and a glycosyltransferase under conditions sufficient to transfer said poly(ethylene glycol) to said EPO peptide, wherein said  
25   glycosyltransferase is a fucosyltransferase.

67. The method of claim 66, wherein said fucosyltransferase is selected from the group consisting of fucosyltransferase I, fucosyltransferase III, fucosyltransferase IV, fucosyltransferase V, fucosyltransferase VI and fucosyltransferase VII.

30           68. A glycoPEGylated EPO peptide made by the method of claim 66.

69. The method of claim 66, wherein said EPO peptide is expressed in a CHO cell.

5 70. A method of treating a mammal having anemia, said method comprising administering to said mammal an EPO peptide having one or more glycans having a glycoconjugate molecule attached to said peptide, wherein said EPO peptide is administered in an amount effective to increase the hematocrit level in said mammal.

10 71. The method of claim 70, wherein said mammal is a human.

72. A method of providing erythropoietin therapy to a mammal, said method comprising administering an effective amount of a glycoPEGylated EPO peptide comprising an EPO peptide and at least one glycan and at least one poly(ethylene glycol) molecule covalently attached to said glycan, wherein said poly(ethylene glycol) molecule is added to  
15 said EPO peptide using a glycosyltransferase, wherein said EPO peptide is administered in an amount effective to increase the hematocrit level in said mammal.

20 73. The method of claim 72, wherein said mammal is a human.

74. A method of treating a mammal having anemia, said method comprising administering to said mammal a glycoPEGylated EPO peptide comprising an EPO peptide and at least one glycan and at least one poly(ethylene glycol) molecule covalently attached to said glycan, wherein said poly(ethylene glycol) molecule is added to said EPO peptide using  
25 a glycosyltransferase, wherein said EPO peptide is administered in an amount effective to increase the hematocrit level in said mammal..

30 75. The method of claim 74, wherein said mammal is a human.




76. The method of claim 75, wherein said anemia is associated with chemotherapy.

77. A method of treating a kidney dialysis patient, said method comprising  
5 administering to said patient a glycoPEGylated EPO peptide comprising an EPO peptide and  
at least one glycan and at least one poly(ethylene glycol) molecule covalently attached to said  
glycan, wherein said poly(ethylene glycol) molecule is added to said EPO peptide using a  
glycosyltransferase, wherein said EPO peptide is administered in an amount effective to  
increase the hematocrit level in said patient.

10

1/498



Trimannosyl core      Trimannosyl core with  
Bisecting GlcNAc

FIG. 1

2/498




FIG. 2

3/498




FIG. 3

4/498




FIG. 4

5/498




FIG. 6

FIG. 5

6/498




FIG. 7

7/498




FIG. 8



8/498




FIG. 9

9/498




FIG. 10

10/498





FIG. 11

11/498



**FIG. 12**

12/498





**Core 1****Core 2****Core 3****Core 4**

FIG. 13



14/498




FIG. 14B

15/498




FIG. 15



16/498




FIG. 16

17/498




FIG. 17

18/498






FIG. 18A

19/498




**FIG. 18B**

20/498



21/498



22/498





FIG. 20

23/498


## O-LINKED OLIGOSACCHARIDES




**FIG. 21** Ser/Thr



24/498



25/498



26/498




FIG. 23C

27/498




FIG. 24

28/498




FIG. 25

29/498




FIG. 26

30/498




FIG. 27

## 31/498

|  |   |
|--|---|
| 12AP1/E5 -- Viventia Biotech               | AI-201 -- AutoImmune                      |
| 1964 -- Aventis                            | AI-301 -- AutoImmune                      |
| 20K growth hormone -- AMUR                 | AIDS vaccine -- ANRS, CIBG, Hesed         |
| 28P6/E6 -- Viventia Biotech                | Biomed, Hollis-Eden, Rome, United         |
| 3-Hydroxyphthaloyl-beta-lactoglobulin --   | Biomedical, American Home Products,       |
| 4-IBB ligand gene therapy --               | Maxygen                                   |
| 64-Cu MAb conjugate TETA-1A3 --            | airway receptor ligand -- IC Innovations  |
| Mallinckrodt Institute of Radiology        | AJvW 2 -- Ajinomoto                       |
| 64-Cu MAb conjugate TETA-cT84.66           | AK 30 NGF -- Alkermes                     |
| 64-Cu Trastuzumab TETA conjugate --        | Albuferon -- Human Genome Sciences        |
| Genentech                                  | albumin -- Biogen, DSM Anti-Infectives,   |
| A 200 -- Amgen                             | Genzyme Transgenics, PPL Therapeutics,    |
| A10255 -- Eli Lilly                        | TranXenoGen, Welfide Corp.                |
| A1PDX -- Hedral Therapeutics               | aldesleukin -- Chiron                     |
| A6 -- Angstrom                             | alefacept -- Biogen                       |
| aaAT-III -- Genzyme                        | Alemtuzumab                               |
| Abciximab -- Centocor                      | Allergy therapy -- ALK-Abello/Maxygen,    |
| ABI.001 -- Atlantic BioPharmaceuticals     | ALK-Abello/RP Scherer                     |
| ABT-828 -- Abbott                          | allergy vaccines -- Allergy Therapeutics  |
| Accutin                                    | Alnidofibatide -- Aventis Pasteur         |
| Actinohivin                                | Alnorine -- SRC VB VECTOR                 |
| activin -- Biotech Australia, Human        | ALP 242 -- Gruenenthal                    |
| Therapeutics, Curis                        | Alpha antitrypsin -- Arriva/Hyland        |
| AD 439 -- Tanox                            | Immuno/ProMetic/Protease Sciences         |
| AD 519 -- Tanox                            | Alpha-1 antitrypsin -- Cutter, Bayer, PPL |
| Adalimumab -- Cambridge Antibody Tech.     | Therapeutics, Profile, ZymoGenetics,      |
| Adenocarcinoma vaccine -- Biomira -- NIS   | Arriva                                    |
| Adenosine deaminase -- Enzond              | Alpha-1 protease inhibitor -- Genzyme     |
| Adenosine A2B receptor antagonists --      | Transgenics, Welfide Corp.                |
| Adenosine Therapeutics                     | Alpha-galactose fusion protein --         |
| ADP-001 -- Axis Genetics                   | Immunomedics                              |
| AF 13948 -- Affymax                        | Alpha-galactosidase A -- Research         |
| Afelimomab -- Knoll                        | Corporation Technologies, Genzyme         |
| AFP-SCAN -- Immunomedics                   | Alpha-glucosidase -- Genzyme, Novazyme    |
| AG 2195 -- Corixa                          | Alpha-lactalbumin                         |
| agalsidase alfa -- Transkaryotic Therapies | Alpha-L-iduronidase -- Transkaryotic      |
| agalsidase beta -- Genzyme                 | Therapies, BioMarin                       |
| AGENT-- Antisoma                           | alteplase -- Genentech                    |
| AI 300 -- AutoImmune                       | alvircept sudotox -- NIH                  |
| AI-101 -- Teva                             | ALX-0600, a GLP-2 agonist -- NPS Allelix  |
| AI-102 -- Teva                             | Corp.                                     |

FIG. 28A



## 32/498

|  |   |
|--|---|
| ALX1-11 --sNPS Pharmaceuticals               | Anti-alphav $\beta$ 3 integrin MAb -- Applied |
| Alzheimer's disease gene therapy             | Molecular Evolution                           |
| AM-133 -- AMRAD                              | Anti-angiogenesis monoclonal antibodies --    |
| Amb a 1 immunostim conj. -- Dynavax          | KS Biomedix/Schering AG                       |
| AMD 3100 -- AnorMED -- NIS                   | Anti-B4 MAb-DC1 conjugate -- ImmunoGen        |
| AMD 3465 -- AnorMED -- NIS                   | Anti-B7 antibody PRIMATIZED -- IDEC           |
| AMD 3465 -- AnorMED -- NIS                   | Anti-B7-1 MAb 16-10A1                         |
| AMD Fab -- Genentech                         | Anti-B7-1 MAb 1G10                            |
| Amediplase -- Menarini, Novartis             | Anti-B7-2 MAb GL-1                            |
| AM-F9  | Anti-B7-2-gelonin immunotoxin --              |
| Amoebiasis vaccine                           | Antibacterials/antifungals --                 |
| Amphiregulin -- Octagene                     | Diversa/IntraBiotics                          |
| anakinra -- Amgen                            | Anti-beta-amyloid monoclonal antibodies --    |
| analgesic -- Nobex                           | Cambridge Antibody Tech., Wyeth-Ayerst        |
| ancestim -- Amgen                            | Anti-BLyS antibodies -- Cambridge             |
| AnergiX.RA -- Corixa, Organon                | Antibody Tech. /Human Genome Sciences         |
| Angiocidin -- InKine                         | Antibody-drug conjugates -- Seattle           |
| angiogenesis inhibitors -- ILEX              | Genetics/Eos                                  |
| AngioMab -- Antisoma                         | Anti-C5 MAb BB5-1 -- Alexion                  |
| Angiopoietins -- Regeneron/Procter &         | Anti-C5 MAb N19-8 -- Alexion                  |
| Gamble                                       | Anti-C8 MAb                                   |
| angiostatin -- EntreMed                      | anticancer cytokines -- BioPulse              |
| Angiostatin/endostatin gene therapy --       | anticancer matrix -- Telios Integra           |
| Genetix Pharmaceuticals                      | Anticancer monoclonal antibodies -- ARIUS,    |
| angiotensin-II, topical -- Maret             | Immunex                                       |
| Anthrax -- EluSys Therapeutics/US Army       | anticancer peptides -- Maxygen, Micrologix    |
| Medical Research Institute                   | Anticancer prodrug Tech. -- Alexion           |
| Anthrax vaccine                              | Antibody Technologies                         |
| Anti platelet-derived growth factor D human  | anticancer Troy-Bodies -- Affite -- Affitech  |
| monoclonal antibodies -- CuraGen             | anticancer vaccine -- NIH                     |
| Anti-17-1A MAb 3622W94 --                    | anticancers -- Epimmune                       |
| GlaxoSmithKline                              | Anti-CCR5/CXCR4 sheep MAb -- KS               |
| Anti-2C4 MAb -- Genentech                    | Biomedix Holdings                             |
| anti-4-1BB monoclonal antibodies -- Bristol- | Anti-CD11a MAb KBA --                         |
| Myers Squibb                                 | Anti-CD11a MAb M17                            |
| Anti-Adhesion Platform Tech. -- Cytovax      | Anti-CD11a MAb TA-3 --                        |
| Anti-adipocyte MAb -- Cambridge Antibody     | Anti-CD11a MAb WT.1 --                        |
| Tech./ObeSys                                 | Anti-CD11b MAb -- Pharmacia                   |
| antiallergics -- Maxygen                     | Anti-CD11b MAb LM2                            |
| antiallergy vaccine -- Acambis               | Anti-CD154 MAb -- Biogen                      |
| Anti-alpha-4-integrin MAb                    | Anti-CD16-anti-CD30 MAb -- Biotest            |

FIG. 28B

## 33/498

|   |   |
|---|---|
| Anti-CD18 MAb -- Pharmacia                                    | Anti-CD4 MAb -- Centocor, IDEC                                      |
| Anti-CD19 MAb B43 --  | Pharmaceuticals, Xenova Group                                       |
| Anti-CD19 MAb -liposomal sodium butyrate conjugate --         | Anti-CD4 MAb 16H5   |
| Anti-CD147  | Anti-CD4 MAb 4162W94 -- GlaxoSmithKline                             |
| Anti-CD19 MAb-saporin conjugate --                            | Anti-CD4 MAb B-F5 -- Diaclone                                       |
| Anti-CD19-dsFv-PE38-immunotoxin --                            | Anti-CD4 MAb GK1-5  |
| Anti-CD2 MAb 12-15 --   | Anti-CD4 MAb KT6  |
| Anti-CD2 MAb B-E2 -- Diaclone                                 | Anti-CD4 MAb OX38   |
| Anti-CD2 MAb OX34 --  | Anti-CD4 MAb PAP conjugate -- Bristol-Myers Squibb                  |
| Anti-CD2 MAb OX54 --  | Anti-CD4 MAb RIB 5-2  |
| Anti-CD2 MAb OX55 --  | Anti-CD4 MAb W3/25  |
| Anti-CD2 MAb RM2-1  | Anti-CD4 MAb YTA 3.1.2  |
| Anti-CD2 MAb RM2-2  | Anti-CD4 MAb YTS 177-9  |
| Anti-CD2 MAb RM2-4  | Anti-CD40 ligand MAb 5c8 -- Biogen                                  |
| Anti-CD20 MAb BCA B20   | Anti-CD40 MAb   |
| Anti-CD20-anti-Fc alpha RI bispecific MAb -- Medarex, Tenovus | Anti-CD40 MAb 5D12 -- Tanox   |
| Anti-CD22 MAb-saporin-6 complex --                            | Anti-CD44 MAb A3D8  |
| Anti-CD3 immunotoxin --                                       | Anti-CD44 MAb GKWA3   |
| Anti-CD3 MAb 145-2C11 -- Pharming                             | Anti-CD44 MAb IM7   |
| Anti-CD3 MAb CD4IgG conjugate -- Genentech                    | Anti-CD44 MAb KM81  |
| Anti-CD3 MAb humanised -- Protein Design, RW Johnson          | Anti-CD44 variant monoclonal antibodies -- Corixa/Hebrew University |
| Anti-CD3 MAb WT32   | Anti-CD45 MAb BC8-I-131   |
| Anti-CD3 MAb-ricin-chain-A conjugate --                       | Anti-CD45RB MAb   |
| Anti-CD3 MAb-xanthine-oxidase conjugate --                    | Anti-CD48 MAb HuLy-m3   |
| Anti-CD30 MAb BerH2 -- Medac                                  | Anti-CD48 MAb WM-63   |
| Anti-CD30 MAb-saporin conjugate                               | Anti-CD5 MAb -- Becton Dickinson                                    |
| Anti-CD30-scFv-ETA'-immunotoxin                               | Anti-CD5 MAb OX19   |
| Anti-CD38 MAb AT13/5  | Anti-CD6 MAb  |
| Anti-CD38 MAb-saporin conjugate                               | Anti-CD7 MAb-PAP conjugate  |
| Anti-CD3-anti-CD19 bispecific MAb                             | Anti-CD7 MAb-ricin-chain-A conjugate                                |
| Anti-CD3-anti-EGFR MAb  | Anti-CD8 MAb -- Amerimmune, Cytodyn, Becton Dickinson               |
| Anti-CD3-anti-interleukin-2-receptor MAb                      | Anti-CD8 MAb 2-43   |
| Anti-CD3-anti-MOV18 MAb -- Centocor                           | Anti-CD8 MAb OX8  |
| Anti-CD3-anti-SCLC bispecific MAb                             | Anti-CD80 MAb P16C10 -- IDEC  |
| Anti-CD4 idiotype vaccine                                     | Anti-CD80 MAb P7C10 -- ID Vaccine                                   |
|   | Anti-CD8-idarubicin conjugate                                       |
|   | Anti-CEA MAb CE-25  |
|   | Anti-CEA MAb MN 14 -- Immunomedics                                  |

FIG. 28C

34/498

|  |  |
|--|--|
| Anti-CEA MAb MN14-PE40 conjugate -- Immunomedics                               | Anti-heparanase human monoclonal antibodies -- Oxford Glycosciences/Medarex          |
| Anti-CEA MAb T84.66-interleukin-2 conjugate                                    | Anti-hepatitis C virus human monoclonal antibodies -- XTL Biopharmaceuticals         |
| Anti-CEA sheep MAb -- KS Biomedix Holdings                                     | Anti-HER-2 antibody gene therapy   |
| Anti-cell surface monoclonal antibodies -- Cambridge Antibody Tech. /Pharmacia | Anti-herpes antibody -- Epicyte  |
| Anti-c-erbB2-anti-CD3 bifunctional MAb -- Otsuka                               | Anti-HIV antibody -- Epicyte   |
| Anti-CMV MAb -- Scotgen  | anti-HIV catalytic antibody -- Hesed Biomed  |
| Anti-complement  | anti-HIV fusion protein -- Idun  |
| Anti-CTLA-4 MAb  | anti-HIV proteins -- Cangene   |
| Anti-EGFR catalytic antibody -- Hesed Biomed                                   | Anti-HM1-24 MAb -- Chugai  |
| anti-EGFR immunotoxin -- IVAX  | Anti-hR3 MAb   |
| Anti-EGFR MAb -- Abgenix   | Anti-Human-Carcinoma-Antigen MAb -- Epicyte  |
| Anti-EGFR MAb 528  | Anti-ICAM-1 MAb -- Boehringer Ingelheim  |
| Anti-EGFR MAb KSB 107 -- KS Biomedix   | Anti-ICAM-1 MAb 1A-29 -- Pharmacia   |
| Anti-EGFR MAb-DM1 conjugate -- ImmunoGen                                       | Anti-ICAM-1 MAb HA58   |
| Anti-EGFR MAb-LA1 --   | Anti-ICAM-1 MAb YN1/1.7.4  |
| Anti-EGFR sheep MAb -- KS Biomedix   | Anti-ICAM-3 MAb ICM3 -- ICOS   |
| Anti-FAP MAb F19-I-131   | Anti-idiotype breast cancer vaccine 11D10  |
| Anti-Fas IgM MAb CH11  | Anti-idiotype breast cancer vaccine ACA14C5 --                                       |
| Anti-Fas MAb Jo2   | Anti-idiotype cancer vaccine -- ImClone Systems/Merck KGaA ImClone, Viventia Biotech |
| Anti-Fas MAb RK-8  | Anti-idiotype cancer vaccine 1A7 -- Titan  |
| Anti-Flt-1 monoclonal antibodies -- ImClone                                    | Anti-idiotype cancer vaccine 3H1 -- Titan  |
| Anti-fungal peptides -- State University of New York                           | Anti-idiotype cancer vaccine TriAb -- Titan  |
| antifungal tripeptides -- BTG  | Anti-idiotype Chlamydia trachomatis vaccine  |
| Anti-ganglioside GD2 antibody-interleukin-2 fusion protein -- Lexigen          | Anti-idiotype colorectal cancer vaccine -- Novartis                                  |
| Anti-GM2 MAb -- Kyowa  | Anti-idiotype colorectal cancer vaccine -- Onyvax                                    |
| Anti-GM-CSF receptor monoclonal antibodies -- AMRAD                            | Anti-idiotype melanoma vaccine -- IDEC Pharmaceuticals                               |
| Anti-gp130 MAb -- Tosoh  | Anti-idiotype ovarian cancer vaccine ACA 125   |
| Anti-HCA monoclonal antibodies -- AltaRex/Epigen                               | Anti-idiotype ovarian cancer vaccine AR54 - AltaRex                                  |
| Anti-hCG antibodies -- Abgenix/AVI BioPharma                                   |  |

**FIG. 28D**

## 35/498

|   |  |
|---|--|
| Anti-idiotypic ovarian cancer vaccine CA-125 – AltaRex, Biomira       | Anti-L-selectin monoclonal antibodies -- Protein Design Labs, Abgenix, Stanford University |
| Anti-IgE catalytic antibody -- Hersed Biomed                          | Anti-MBL monoclonal antibodies -- Alexion/Brigham and Women's Hospital                     |
| Anti-IgE MAb E26 -- Genentech   | Anti-MHC monoclonal antibodies   |
| Anti-IGF-1 MAb  | Anti-MIF antibody humanised – IDEC, Cytokine PharmaSciences                                |
| anti-inflammatory -- GeneMax  | Anti-MRSA/VRSA sheep MAb -- KS Biomedix Holdings   |
| anti-inflammatory peptide -- BTG                                      | Anti-mu MAb -- Novartis  |
| anti-integrin peptides -- Burnha                                      | Anti-MUC-1 MAb   |
| Anti-interferon-alpha-receptor MAb 64G12 -- Pharma Pacific Management | Anti-MUC 18  |
| Anti-interferon-gamma MAb -- Protein Design Labs                      | Anti-Nogo-A MAb IN1  |
| Anti-interferon-gamma polyclonal antibody - Advanced Biotherapy       | Anti-nuclear autoantibodies -- Procyon   |
| Anti-interleukin-10 MAb --  | Anti-ovarian cancer monoclonal antibodies - Dompe  |
| Anti-interleukin-12 MAb --  | Anti-p185 monoclonal antibodies  |
| Anti-interleukin-1-beta polyclonal antibody -- R&D Systems            | Anti-p43 MAb   |
| Anti-interleukin-2 receptor MAb 2A3                                   | Antiparasitic vaccines   |
| Anti-interleukin-2 receptor MAb 33B3-1 -- Immunotech                  | Anti-PDGF/bFGF sheep MAb -- KS Biomedix  |
| Anti-interleukin-2 receptor MAb ART-18                                | Anti-properdin monoclonal antibodies -- Abgenix/Gliatech                                   |
| Anti-interleukin-2 receptor MAb LO-Tact-1                             | Anti-PSMA (prostate specific membrane antigen)   |
| Anti-interleukin-2 receptor MAb Mikbeta1                              | Anti-PSMA MAb J591 -- BZL Biologics  |
| Anti-interleukin-2 receptor MAb NDS61                                 | Anti-Rev MAb gene therapy --   |
| Anti-interleukin-4 MAb 11B11  | Anti-RSV antibodies -- Epicyte, Intracell  |
| Anti-interleukin-5 MAb -- Wallace Laboratories                        | Anti-RSV monoclonal antibodies -- Medarex/MedImmune, Applied Molecular Evolution/MedImmune |
| Anti-interleukin-6 MAb -- Centocor, Diaclone, Pharmadigm              | Anti-RSV MAb, inhalation -- Alkermes/MedImmune   |
| Anti-interleukin-8 MAb -- Abgenix                                     | Anti-RT gene therapy   |
| Anti-interleukin-8 MAb -- Xenotech                                    | Antisense K-ras RNA gene therapy   |
| Anti-JL1 MAb  | Anti-SF-25 MAb   |
| Anti-Klebsiella sheep MAb -- KS Biomedix Holdings                     | Anti-sperm antibody -- Epicyte   |
| Anti-Laminin receptor MAb-liposomal doxorubicin conjugate             | Anti-Tac(Fv)-PE38 conjugate  |
| Anti-LCG MAb -- Cytoclonal  | Anti-TAPA/CD81 MAb AMP1  |
| Anti-lipopolysaccharide MAb -- VitaResc                               | Anti-tat gene therapy  |

FIG. 28E

## 36/498

|  |   |
|--|---|
| Anti-TCR-alphabeta MAb H57-597   | AOP-RANTES -- Senetek   |
| Anti-TCR-alphabeta MAb R73   | Apan-CH -- Praecis Pharmaceuticals  |
| Anti-tenascin MAb BC-4-I-131   | APC-8024 -- Demegen   |
| Anti-TGF-beta human monoclonal antibodies -- Cambridge Antibody Tech., Genzyme     | ApoA-1 -- Milano, Pharmacia   |
| Anti-TGF-beta MAb 2G7 -- Genentech   | Apogen -- Alexion   |
| Antithrombin III -- Genzyme Transgenics, Aventis, Bayer, Behringwerke, CSL, Myriad | apolipoprotein A1 -- Avanir   |
| Anti-Thy1 MAb  | Apolipoprotein E -- Bio-Tech. General   |
| Anti-Thy1.1 MAb  | Applaggin -- Biogen   |
| Anti-tissue factor/factor VIIA sheep MAb -- KS Biomedix                            | aprotinin -- ProdiGene  |
| Anti-TNF monoclonal antibodies -- Centocor, Chiron, Peptech, Pharacia, Serono      | APT-070C -- AdProTech   |
| Anti-TNF sheep MAb -- KS Biomedix Holdings   | AR 177 -- Aronex Pharmaceuticals  |
| Anti-TNFalpha MAb -- Genzyme   | AR 209 -- Aronex Pharmaceuticals, Antigenics                                  |
| Anti-TNFalpha MAb B-C7 -- Diaclone   | AR545C  |
| Anti-tooth decay MAb -- Planet BioTech.  | ARGENT gene delivery systems -- ARIAD   |
| Anti-TRAIL receptor-1 MAb -- Takeda  | Arresten  |
| Antitumour RNases -- NIH   | ART-123 -- Asahi Kasei  |
| Anti-VCAM MAb 2A2 -- Alexion   | arylsulfatase B -- BioMarin   |
| Anti-VCAM MAb 3F4 -- Alexion   | Arylsulfatase B, Recombinant human -- BioMarin                                |
| Anti-VCAM-1 MAb  | AS 1051 -- Ajinomoto  |
| Anti-VEC MAb -- ImClone  | ASI-BCL -- Intracell  |
| Anti-VEGF MAb -- Genentech   | Asparaginase - Merck  |
| Anti-VEGF MAb 2C3  | ATL-101 -- Alizyme  |
| Anti-VEGF sheep MAb -- KS Biomedix Holdings  | Atrial natriuretic peptide -- Pharis  |
| Anti-VLA-4 MAb HP1/2 -- Biogen   | Aurintricarboxylic acid-high molecular weight                                 |
| Anti-VLA-4 MAb PS/2  | Autoimmune disorders -- GPC Biotech/MorphoSys                                 |
| Anti-VLA-4 MAb R1-2  | Autoimmune disorders and transplant rejection -- Bristol-Myers Squibb/Genzyme |
| Anti-VLA-4 MAb TA-2  | Tra   |
| Anti-VAP-1 human MAb   | Autoimmune disorders/cancer -- Abgenix/Chiron, CuraGen                        |
| Anti-VRE sheep MAb -- KS Biomedix Holdings   | Autotaxin   |
| ANUP -- TranXenoGen  | Avicidin -- NeoRx   |
| ANUP-1 -- Pharis   | axogenesis factor-1 -- Boston Life Sciences                                   |
|  | Axokine -- Regeneron  |
|  | B cell lymphoma vaccine -- Biomira  |
|  | B7-1 gene therapy --  |
|  | BABS proteins -- Chiron   |

FIG. 28F

37/498

|  |   |
|--|---|
| BAM-002 -- Novelos Therapeutics              | BMP 2 -- Genetics Institute/Medtronic-      |
| Basiliximab (anti CD25 MAb) -- Novartis      | Sofamor Danek, Genetics Institute/          |
| Bay-16-9996 -- Bayer                         | Collagenesis, Genetics                      |
| Bay-39-9437 -- Bayer                         | Institute/Yamanouch                         |
| Bay-50-4798 -- Bayer                         | BMP 2 gene therapy                          |
| BB-10153 -- British Biotech                  | BMP 52 -- Aventis Pasteur, Biopharm         |
| BBT-001 -- Bolder BioTech.                   | BMP-2 -- Genetics Institute                 |
| BBT-002 -- Bolder BioTech.                   | BMS 182248 -- Bristol-Myers Squibb          |
| BBT-003 -- Bolder BioTech.                   | BMS 202448 -- Bristol-Myers Squibb          |
| BBT-004 -- Bolder BioTech.                   | bone growth factors -- IsoTis               |
| BBT-005 -- Bolder BioTech.                   | BPC-15 -- Pfizer                            |
| BBT-006 -- Bolder BioTech.                   | brain natriuretic peptide --                |
| BBT-007 -- Bolder BioTech.                   | Breast cancer -- Oxford                     |
| BCH-2763 -- Shire                            | GlycoSciences/Medarex                       |
| BCSF -- Millenium Biologix                   | Breast cancer vaccine -- Therion Biologics, |
| BDNF -- Regeneron -- Amgen                   | Oregon                                      |
| Becaplermin -- Johnson & Johnson, Chiron     | BSSL -- PPL Therapeutics                    |
| Bectumomab -- Immunomedics                   | BST-2001 -- BioStratum                      |
| Beriplast -- Aventis                         | BST-3002 -- BioStratum                      |
| Beta-adrenergic receptor gene therapy --     | BTI 322 --                                  |
| University of Arkansas                       | butyrylcholinesterase -- Shire              |
| bFGF -- Scios                                | C 6822 -- COR Therapeutics                  |
| BI 51013 -- Behringwerke AG                  | C1 esterase inhibitor -- Pharming           |
| BIBH 1 -- Boehringer Ingelheim               | C3d adjuvant -- AdProTech                   |
| BIM-23190 -- Beaufour-Ipsen                  | CAB-2.1 -- Millennium                       |
| birch pollen immunotherapy -- Pharmacia      | calcitonin -- Inhale Therapeutics Systems,  |
| bispecific fusion proteins -- NIH            | Aventis, Genetronics, TranXenoGen,          |
| Bispecific MAb 2B1 -- Chiron                 | Unigene, Rhone Poulenc Rohrer               |
| Bitistatin                                   | calcitonin -- oral -- Nobex, Emisphere,     |
| BIWA 4 -- Boehringer Ingelheim               | Pharmaceutical Discovery                    |
| blood substitute -- Northfield, Baxter Intl. | Calcitonin gene-related peptide -- Asahi    |
| BLP-25 -- Biomira                            | Kasei -- Unigene                            |
| BLS-0597 -- Boston Life Sciences             | calcitonin, human -- Suntory                |
| BLyS -- Human Genome Sciences                | calcitonin, nasal -- Novartis, Unigene      |
| BLyS radiolabelled -- Human Genome           | calcitonin, Panoderm -- Elan                |
| Sciences                                     | calcitonin, Peptitrol -- Shire              |
| BM 06021 -- Boehringer Mannheim              | calcitonin, salmon -- Therapicon            |
| BM-202 -- BioMarin                           | calin -- Biopharm                           |
| BM-301 -- BioMarin                           | Calphobindin I                              |
| BM-301 -- BioMarin                           | calphobindin I -- Kowa                      |
| BM-302 -- BioMarin                           | calreticulin -- NYU                         |

**FIG. 28G**

## 38/498

|  |   |
|--|---|
| Campath-1G                                 | CD4 fusion toxin -- Senetek             |
| Campath-1M                                 | CD4 IgG -- Genentech                    |
| cancer therapy -- Cangene                  | CD4 receptor antagonists --             |
| cancer vaccine -- Aixlie, Aventis Pasteur, | Pharmacoepia/Progenics                  |
| Center of Molecular Immunology, YM         | CD4 soluble -- Progenics                |
| BioSciences, Cytos, Genzyme,               | CD4, soluble -- Genzyme Transgenics     |
| Transgenics, GlobelImmune, Igeneon,        | CD40 ligand -- Immunex                  |
| ImClone, Virogenetics, InterCell, Iomai,   | CD4-ricin chain A -- Genentech          |
| Jenner Biotherapies, Memorial Sloan-       | CD59 gene therapy -- Alexion            |
| Kettering Cancer Center, Sydney Kimmel     | CD8 TIL cell therapy -- Aventis Pasteur |
| Cancer Center, Novavax, Protein            | CD8, soluble -- Avidex                  |
| Sciences, Argonex, SIGA                    | CD95 ligand -- Roche                    |
| Cancer vaccine ALVAC-CEA B7.1 --           | CDP 571 -- Celltech                     |
| Aventis Pasteur/Therion Biologics          | CDP 850 -- Celltech                     |
| Cancer vaccine CEA-TRICOM -- Aventis       | CDP-860 (PEG-PDGF MAb) -- Celltech      |
| Pasteur/Therion Biologics                  | CDP 870 -- Celltech                     |
| Cancer vaccine gene therapy -- Cantab      | CDS-1 -- Ernest Orlando                 |
| Pharmaceuticals                            | Cedelizumab -- Ortho-McNeil             |
| Cancer vaccine HER-2/neu -- Corixa         | Cetermin -- Insmad                      |
| Cancer vaccine THERATOPE -- Biomira        | CETP vaccine -- Avant                   |
| cancer vaccine, PolyMASC -- Valentis       | Cetrorelix                              |
| Candida vaccine -- Corixa, Inhibitex       | Cetuximab                               |
| Canstatin -- ILEX                          | CGH 400 -- Novartis                     |
| CAP-18 -- Panorama                         | CGP 42934 -- Novartis                   |
| Cardiovascular gene therapy -- Collateral  | CGP 51901 -- Tanox                      |
| Therapeutics                               | CGRP -- Unigene                         |
| carperitide -- Suntory                     | CGS 27913 -- Novartis                   |
| Casocidin-1 -- Pharis                      | CGS 32359 -- Novartis                   |
| CAT 152 -- Cambridge Antibody Tech.        | Chagas disease vaccine -- Corixa        |
| CAT 192 -- Cambridge Antibody Tech.        | chemokines -- Immune Response           |
| CAT 213 -- Cambridge Antibody Tech.        | CHH 380 -- Novartis                     |
| Catalase-- Enzon                           | chitinase -- Genzyme, ICOS              |
| Cat-PAD -- Circassia                       | Chlamydia pneumoniae vaccine -- Antex   |
| CB 0006 -- Celltech                        | Biologics                               |
| CCK(27-32)-- Akzo Nobel                    | Chlamydia trachomatis vaccine -- Antex  |
| CCR2-64I -- NIH                            | Biologics                               |
| CD, Procept -- Paligent                    | Chlamydia vaccine -- GlaxoSmithKline    |
| CD154 gene therapy                         | Cholera vaccine CVD 103-HgR -- Swiss    |
| CD39 -- Immunex                            | Serum and Vaccine Institute Berne       |
| CD39-L2 -- Hyseq                           | Cholera vaccine CVD 112 -- Swiss Serum  |
| CD39-L4 -- Hyseq                           | and Vaccine Institute Berne             |

FIG. 28H

## 39/498

|  |  |
|--|--|
| Cholera vaccine inactivated oral -- SBL Vaccin                                   | CRL 1605 -- CytRx  |
| Chrysalin -- Chrysalis BioTech.  | CS-560 -- Sankyo   |
| CI-782 -- Hitachi Kase   | CSF -- ZymoGenetics  |
| Ciliary neurotrophic factor -- Fidia, Roche                                      | CSF-G -- Hangzhou, Dong-A, Hanmi   |
| CIM project -- Active Biotech  | CSF-GM -- Cangene, Hunan, LG Chem  |
| CL 329753 -- Wyeth-Ayerst  | CSF-M -- Zarix   |
| CL22, Cobra -- ML Laboratories   | CT 1579 -- Merck Frosst  |
| Clenoliximab -- IDEC   | CT 1786 -- Merck Frosst  |
| Clostridium difficile antibodies -- Epicyte                                      | CT-112 <sup>^</sup> -- BTG   |
| clotting factors -- Octagene   | CTB-134L -- Xenova   |
| CMB 401 -- Celltech  | CTC-111 -- Kaketsuken  |
| CNTF -- Sigma-Tau  | CTGF -- FibroGen   |
| Cocaine abuse vaccine -- Cantab, ImmuLogic, Scripps                              | CTLA4-Ig -- Bristol-Myers Squibb   |
| coccidiomycosis vaccine -- Arizo   | CTLA4-Ig gene therapy --   |
| collagen -- Type I -- Pharming   | CTP-37 -- AVI BioPharma  |
| Collagen formation inhibitors -- FibroGen  | C-type natriuretic peptide -- Suntory  |
| Collagen/hydroxyapatite/bone growth factor -- Aventis Pasteur, Biopharm, Orquest | CVS 995 -- Corvas Intl.  |
| collagenase -- BioSpecifics  | CX 397 -- Nikko Kyodo  |
| Colorectal cancer vaccine -- Wistar Institute                                    | CY 1747 -- Epimmune  |
| Component B, Recombinant -- Serono   | CY 1748 -- Epimmune  |
| Connective tissue growth factor inhibitors -- FibroGen/Taisho                    | Cyanovirin-N   |
| Contortrostatin  | Cystic fibrosis therapy -- CBR/IVAX  |
| contraceptive vaccine -- Zonagen   | CYT 351  |
| Contraceptive vaccine hCG  | cytokine Traps -- Regeneron  |
| Contraceptive vaccine male reversible -- IMMUCON                                 | cytokines -- Enzon, Cytoclonal   |
| Contraceptive vaccine zona pellucida -- Zonagen                                  | Cytomegalovirus glycoprotein vaccine -- Chiron, Aquila Biopharmaceuticals, Aventis Pasteur, Virogenetics |
| Copper-64 labelled MAb TETA-1A3 -- NCI   | Cytomegalovirus vaccine live -- Aventis Pasteur  |
| Coralyne   | Cytosine deaminase gene therapy -- GlaxoSmithKline   |
| Corsevin M   | DA-3003 -- Dong-A  |
| C-peptide analogues -- Schwarz   | DAB389interleukin-6 -- Senetek   |
| CPI-1500 -- Consensus  | DAB389interleukin-7  |
| CRF -- Neurobiological Tech.   | DAC:GLP-2 -- ConjuChem, Inc.   |
| cRGDfV pentapeptide --   | Daclizumab (anti-IL2R MAb) -- Protein Design Labs  |
| CRL 1095 -- CytRx  | DAMP <sup>^</sup> -- Incyte Genomics   |
| CRL 1336 -- CytRx  | Daniplestim -- Pharmacia   |
|  | darbepoetin alfa -- Amgen  |

FIG. 28I



## 40/498

|   |  |
|---|--|
| DBI-3019 -- Diabetogen  | dural graft matrix -- Integra                      |
| DCC -- Genzyme  | Duteplase -- Baxter Intl.                          |
| DDF -- Hyseq  | DWP-401 -- Daewoong                                |
| decorin -- Integra, Telios  | DWP-404 -- Daewoong                                |
| defensins -- Large Scale Biology                                    | DWP-408 -- Daewoong                                |
| DEGR-VIIa   | Dx 88 (Epi-KAL2) -- Dyax                           |
| Delimmunised antibody 3B6/22 AGEN                                   | Dx 890 (elastin inhibitors) -- Dyax                |
| Deimmunised anti-cancer antibodies -- Biovation/Viragen             | E coli O157 vaccine -- NIH                         |
| Dendroamide A   | E21-R -- BresaGen                                  |
| Dengue vaccine -- Bavarian Nordic, Merck                            | Eastern equine encephalitis virus vaccine --       |
| denileukin diftitox -- Ligand                                       | Echicetin --                                       |
| DES-1101 -- Desmos  | Echinhibin 1 --                                    |
| desirudin -- Novartis   | Echistatin -- Merck                                |
| desmopressin -- Unigene   | Echitamine --                                      |
| Desmoteplase -- Merck, Schering AG                                  | Ecromeximab -- Kyowa Hakko                         |
| Destabilase   | EC-SOD -- PPL Therapeutics                         |
| Diabetes gene therapy -- DeveloGen, Pfizer                          | Eculizumab (5G1.1) -- Alexion                      |
| Diabetes therapy -- Crucell   | EDF -- Ajinomoto                                   |
| Diabetes type 1 vaccine -- Diamyd Therapeutics                      | EDN derivative -- NIH                              |
| DiaCIM -- YM BioSciences  | EDNA -- NIH  |
| dialytic oligopeptides -- Research Corp                             | Edobacomab -- XOMA                                 |
| Diamyd -- Diamyd Therapeutics                                       | Edrecolomab -- Centocor                            |
| DiaPep227 -- Pepgen   | EF 5077  |
| DiavaX -- Corixa  | Efalizumab -- Genentech                            |
| Digoxin MAb -- Glaxo  | EGF fusion toxin -- Seragen, Ligand                |
| Diphtheria tetanus pertussis-hepatitis B vaccine -- GlaxoSmithKline | EGF-P64k vaccine -- Center of Molecular Immunology |
| DIR therapy -- Solis Therapeutics --                                | EL 246 -- LigoCyte                                 |
| DNase -- Genentech  | elastase inhibitor -- Synergen                     |
| Dornase alfa -- Genentech   | elcatonin -- Therapicon                            |
| Dornase alfa, inhalation -- Genentech                               | EMD 72000 -- Merck KGaA                            |
| Doxorubicin-anti-CEA MAb conjugate -- Immunomedics                  | Emdogain -- BIORA                                  |
| DP-107 -- Trimeris  | emfilermin -- AMRAD                                |
| drotrecogin alfa -- Eli Lilly                                       | Emoctakin -- Novartis                              |
| DTctGMCSF   | enamel matrix protein -- BIORA                     |
| DTP-polio vaccine -- Aventis Pasteur                                | Endo III -- NYU                                    |
| DU 257-KM231 antibody conjugate -- Kyowa                            | endostatin -- EntreMed, Pharis                     |
|   | Enhancins -- Micrologix                            |
|   | Enlimomab -- Isis Pharm.                           |
|   | Enoxaparin sodium -- Pharmuka                      |

FIG. 28J

## 41/498

|   |   |
|---|---|
| enzyme linked antibody nutrient depletion therapy -- KS Biomedix Holdings   | Factor IX gene therapy -- Cell Genesys  |
| Eosinophil-derived neutralizing agent -- EP-51216 -- Asta Medica  | Factor VII -- Novo Nordisk, Bayer, Baxter Intl.   |
| EP-51389 -- Asta Medica   | Factor VIIa -- PPL Therapeutics, ZymoGenetics   |
| EPH family ligands -- Regeneron   | Factor VIII -- Bayer Genentech, Beaufour-Ipsen, CLB, Inex, Octagen, Pharmacia, Pharming |
| Epidermal growth factor -- Hitachi Kasei, Johnson & Johnson   | Factor VIII -- PEGylated -- Bayer   |
| Epidermal growth factor fusion toxin -- Senetek   | Factor VIII fragments -- Pharmacia  |
| Epidermal growth factor-genistein -- EPI-HNE-4 -- Dyax  | Factor VIII gene therapy -- Targeted Genetics   |
| EPI-KAL2 -- Dyax  | Factor VIII sucrose formulation -- Bayer, Genentech                                     |
| Epoetin-alfa -- Amgen, Dragon Pharmaceuticals, Nanjing Huaxin   | Factor VIII-2 -- Bayer  |
| Epratuzumab -- Immunomedics   | Factor VIII-3 -- Bayer  |
| Epstein-Barr virus vaccine -- Aviron/SmithKline Beecham, Bioresearch  | Factor Xa inhibitors -- Merck, Novo Nordisk, Mochida                                    |
| Eptacog alfa -- Novo Nordisk  | Factor XIII -- ZymoGenetics   |
| Eptifibatide -- COR Therapeutics  | Factors VIII and IX gene therapy -- Genetics Institute/Targeted Genetics                |
| erb-38 --   | Famoxin -- Genset   |
| Erlizumab -- Genentech  | Fas (delta) TM protein -- LXR BioTech.  |
| erythropoietin -- Alkermes, ProLease, Dong-A, Elanex, Genetics Institute, LG Chem, Protein Sciences, Serono, Snow Brand, SRC VB VECTOR, Transkaryotic Therapies | Fas TR -- Human Genome Sciences   |
| Erythropoietin Beta -- Hoffman La Roche   | Felvizumab -- Scotgen   |
| Erythropoietin/Epoetin alfa -- Chugai   | FFR-VIIa -- Novo Nordisk  |
| Escherichia coli vaccine -- North American Vaccine, SBL Vaccin, Swiss Serum and Vaccine Institute Berne   | FG-001 -- F-Gene  |
| etanercept -- Immunex   | FG-002 -- F-Gene  |
| examorelin -- Mediolanum  | FG-004 -- F-Gene  |
| Exendin 4 -- Amylin   | FG-005 -- F-Gene  |
| exonuclease VII   | FGF + fibrin -- Repair  |
| F 105 -- Centocor   | Fibrimage -- Bio-Tech. General  |
| F-992 -- Fornix   | fibrin-binding peptides -- ISIS Innovation  |
| Factor IX -- Alpha Therapeutics, Welfide Corp., CSL, enetics Institute/AHP, Pharmacia, PPL Therapeutics   | fibrinogen -- PPL Therapeutics, Pharming  |
|   | fibroblast growth factor -- Chiron, NYU, Ramot, ZymoGenetics                            |
|   | fibrolase conjugate -- Schering AG  |
|   | Filgrastim -- Amgen   |
|   | filgrastim -- PDA modified -- Xencor  |
|   | FLT-3 ligand -- Immunex   |
|   | FN18 CRM9 --  |

**FIG. 28K**

## 42/498

|  |   |
|--|---|
| follistatin -- Biotech Australia, Human Therapeutics                     | Glucocerebrosidase -- Genzyme   |
| follitropin alfa -- Alkermes, ProLease, PowderJect, Serono, Akzo Nobel   | glutamate decarboxylase -- Genzyme Transgenics                            |
| Follitropin Beta -- Bayer, Organon                                       | Glycoprotein S3 -- Kureha   |
| FP 59  | GM-CSF -- Immuhex   |
| FSH -- Ferring   | GM-CSF tumour vaccine -- PowderJect                                       |
| FSH + LH -- Ferring  | GnRH immunotherapeutic -- Protherics                                      |
| F-spondin -- CeNeS   | Goserelin (LhRH antagonist) -- AstraZeneca                                |
| fusion protein delivery system -- UAB Research Foundation                | gp75 antigen -- ImClone   |
| fusion toxins -- Boston Life Sciences                                    | gp96 -- Antigenics  |
| G 5598 -- Genentech  | GPI 0100 -- Galenica  |
| GA-II -- Transkaryotic Therapies   | GR 4991W93 -- GlaxoSmithKline   |
| Gamma-interferon analogues -- SRC VB VECTOR                              | Granulocyte colony-stimulating factor -- Dong-A                           |
| Ganirelix -- Roche   | Granulocyte colony-stimulating factor conjugate                           |
| gastric lipase -- Meristem   | grass allergy therapy -- Dynavax  |
| Gavilimomab --   | GRF1-44 -- ICN  |
| G-CSF -- Amgen, SRC VB VECTOR  | Growth Factor -- Chiron, Atrigel, Atrix, Innogenetics, ZymoGenetics, Novo |
| GDF-1 -- CeNeS   | growth factor peptides -- Biotherapeutics                                 |
| GDF-5 -- Biopharm  | growth hormone -- LG Chem   |
| GDNF (glial derived neurotrophic factor) -- Amgen                        | growth hormone, Recombinant human -- Serono                               |
| gelsolin -- Biogen   | GT 4086 -- Gliatech   |
| Gemtuzumab ozogamicin -- Celltech  | GW 353430 -- GlaxoSmithKline  |
| Gene-activated epoetin-alfa -- Aventis Pharma -- Transkaryotic Therapies | GW-278884 -- GlaxoSmithKline  |
| Glanzmann thrombasthenia gene therapy --                                 | H 11 -- Viventia Biotech  |
| Glatiramer acetate -- Yeda   | H5N1 influenza A virus vaccine -- Protein Sciences                        |
| glial growth factor 2 -- CeNeS   | haemoglobin -- Biopure  |
| GLP-1 -- Amylin, Suntory, TheraTech, Watson                              | haemoglobin 3011, Recombinant -- Baxter Healthcare                        |
| GLP-1 peptide analogues -- Zealand Pharmaceuticals                       | haemoglobin crosfumaril -- Baxter Intl.                                   |
| GLP-2 -- Novo Nordisk, Ontario, Inc., Suntory Limited                    | haemoglobin stabilized -- Ajinomoto                                       |
| glucagon -- Eli Lilly, ZymoGenetics                                      | haemoglobin, recombinant -- Apex  |
| Glucagon-like peptide-1 7-36 amide -- Suntory                            | HAF -- Immune Response  |
| Glucogen-like peptide -- Amylin  | Hantavirus vaccine  |
|  | HB 19   |
|  | HBNF -- Regeneron   |
|  | HCC-1 -- Pharis   |

FIG. 28L

## 43/498

|  |  |
|--|--|
| hCG -- Milkhaus                              | Herpes simplex glycoprotein DNA vaccine -- |
| hCG vaccine -- Zonagen                       | Merck, Wyeth-Lederle Vaccines-Malvern,     |
| HE-317 -- Hollis-Eden Pharmaceuticals        | Genentech, GlaxoSmithKline, Chiron,        |
| Heat shock protein cancer and influenza      | Takeda                                     |
| vaccines -- StressGen                        | Herpes simplex vaccine -- Cantab           |
| Helicobacter pylori vaccine -- Acambis,      | Pharmaceuticals, CEL-SCI, Henderson        |
| AstraZeneca/CSL, Chiron, Provalis            | Morley                                     |
| Helistat-G -- GalaGen                        | Herpes simplex vaccine live -- ImClone     |
| Hemolink -- Hemosol                          | Systems/Wyeth-Lederle, Aventis Pasteur     |
| hepapoietin -- Snow Brand                    | HGF derivatives -- Dompe                   |
| heparanase -- InSight                        | hIAPP vaccine -- Crucell                   |
| heparinase I -- Ibex                         | Hib-hepatitis B vaccine -- Aventis Pasteur |
| heparinase III -- Ibex                       | HIC 1                                      |
| Hepatitis A vaccine -- American Biogenetic   | HIP-- Altachem                             |
| Sciences                                     | Hirudins -- Biopharma, Cangene, Dongkook,  |
| Hepatitis A vaccine inactivated              | Japan Energy Corporation, Pharmacia        |
| Hepatitis A vaccine Nothav -- Chiron         | Corporation, SIR International, Sanofi-    |
| Hepatitis A-hepatitis B vaccine --           | Synthelabo, Sotragene, Rhein Biotech       |
| GlaxoSmithKline                              | HIV edible vaccine -- ProdiGene            |
| hepatitis B therapy -- Tripep                | HIV gp120 vaccine -- Chiron, Ajinomoto,    |
| Hepatitis B vaccine -- Amgen, Chiron SpA,    | GlaxoSmithKline, ID Vaccine, Progenics,    |
| Meiji Milk, NIS, Prodeva, PowderJect,        | VaxGen                                     |
| Rhein Biotech                                | HIV gp120 vaccine gene therapy --          |
| Hepatitis B vaccine recombinant -- Evans     | HIV gp160 DNA vaccine -- PowderJect,       |
| Vaccines, Epitec Combiotech, Genentech,      | Aventis Pasteur, Oncogen, Hyland           |
| MedImmune, Merck Sharp & Dohme,              | Immuno, Protein Sciences                   |
| Rhein Biotech, Shantha Biotechnics,          | HIV gp41 vaccine -- Panacos                |
| Vector, Yeda                                 | HIV HGP-30W vaccine -- CEL-SCI             |
| Hepatitis B vaccine recombinant TGP 943 --   | HIV immune globulin -- Abbott, Chiron      |
| Takeda                                       | HIV peptides -- American Home Products     |
| Hepatitis C vaccine -- Bavarian Nordic,      | HIV vaccine -- Applied bioTech., Axis      |
| Chiron, Innogenetics Acambis,                | Genetics, Biogen, Bristol-Myers Squibb,    |
| Hepatitis D vaccine -- Chiron Vaccines       | Genentech, Korea Green Cross, NIS,         |
| Hepatitis E vaccine recombinant --           | Oncogen, Protein Sciences Corporation,     |
| Genelabs/GlaxoSmithKline, Novavax            | Terumo, Tonen Corporation, Wyeth-          |
| hepatocyte growth factor -- Panorama,        | Ayerst, Wyeth-Lederle Vaccines-Malvern,    |
| Sosei  | Advanced BioScience Laboratories,          |
| hepatocyte growth factor kringle fragments - | Bavarian Nordic, Bavarian Nordic/Statens   |
| - EntreMed                                   | Serum Institute, GeneCure, Immune          |
| Her-2/Neu peptides -- Corixa                 | Response, Progenics, Therion Biologics,    |
|  | United Biomedical, Chiron                  |

FIG. 28M

44/498

|  |  |
|--|--|
| HIV vaccine vCP1433 -- Aventis Pasteur   | Human monoclonal antibodies --             |
| HIV vaccine vCP1452 -- Aventis Pasteur   | Medarex/Northwest Biotherapeutics,         |
| HIV vaccine vCP205 -- Aventis Pasteur    | Medarex/Seattle Genetics                   |
| HL-9 -- American BioScience              | human netrin-1 -- Exelixis                 |
| HM-9239 -- Cytran                        | human papillomavirus antibodies -- Epicyte |
| HML-103 -- Hemosol                       | Human papillomavirus vaccine -- Biotech    |
| HML-104 -- Hemosol                       | Australia, IDEC, StressGen                 |
| HML-105 -- Hemosol                       | Human papillomavirus vaccine MEDI 501 --   |
| HML-109 -- Hemosol                       | MedImmune/GlaxoSmithKline                  |
| HML-110 -- Hemosol                       | Human papillomavirus vaccine MEDI          |
| HML-121 -- Hemosol                       | 503/MEDI 504 --                            |
| hNLP -- Pharis                           | MedImmune/GlaxoSmithKline                  |
| Hookworm vaccine                         | Human papillomavirus vaccine TA-CIN --     |
| host-vector vaccines -- Henogen          | Cantab Pharmaceuticals                     |
| HPM 1 -- Chugai                          | Human papillomavirus vaccine TA-HPV --     |
| HPV vaccine -- MediGene                  | Cantab Pharmaceuticals                     |
| HSA -- Meristem                          | Human papillomavirus vaccine TH-GW --      |
| HSF -- StressGen                         | Cantab/GlaxoSmithKline                     |
| HSP carriers --Weizmann, Yeda, Peptor    | human polyclonal antibodies -- Biosite/Eos |
| HSPPC-70 -- Antigenics                   | BioTech./ Medarex                          |
| HSPPC-96, pathogen-derived -- Antigenics | human type II anti factor VIII monoclonal  |
| HSV 863 -- Novartis                      | antibodies -- ThromboGenics                |
| HTLV-I DNA vaccine                       | humanised anti glycoprotein Ib murine      |
| HTLV-I vaccine                           | monoclonal antibodies -- ThromboGenics     |
| HTLV-II vaccine -- Access                | HumaRAD -- Intracell                       |
| HU 901 -- Tanox                          | HuMax EGFR -- Genmab                       |
| Hu23F2G -- ICOS                          | HuMax-CD4 -- Medarex                       |
| HuHMFG1                                  | HuMax-IL15 -- Genmab                       |
| HumaLYM -- Intracell                     | HYB 190 -- Hybridon                        |
| Human krebs statika -- Yamanouchi        | HYB 676 -- Hybridon                        |
| human monoclonal antibodies --           | I-125 MAb A33 -- Celltech                  |
| Abgenix/Biogen, Abgenix/ Corixa,         | Ibritumomab tiuxetan -- IDEC               |
| Abgenix/Immunex, Abgenix/Lexicon,        | IBT-9401 -- Ibex                           |
| Abgenix/ Pfizer, Athersys/Medarex,       | IBT-9402 -- Ibex                           |
| Biogen/MorphoSys, CAT/Searle,            | IC 14 -- ICOS                              |
| Centocor/Medarex, Corixa/Kirin Brewery,  | Idarubicin anti-Ly-2.1 --                  |
| Corixa/Medarex, Eos BioTech./Medarex,    | IDEC 114 -- IDEC                           |
| Eos/Xenerex, Exelixis/Protein Design     | IDEC 131 -- IDEC                           |
| Labs, ImmunoGen/ Raven, Medarex/         | IDEC 152 -- IDEC                           |
| B.Twelve, MorphoSys/ImmunoGen, XTL       | IDM 1 -- IDM                               |
| Biopharmaceuticals/Dyax,                 | IDPS -- Hollis-Eden Pharmaceuticals        |

**FIG. 28N**

## 45/498

|  |   |
|--|---|
| iduronate-2-sulfatase -- Transkaryotic Therapies             | insulin -- AutoImmune, Altea, Biobras, BioSante, Bio-Tech. General, Chong Kun Dang, Emisphere, Flamel, Provalis, Rhein Biotech, TranXenoGen |
| IGF/IBP-2-13 -- Pharis                                       | insulin (bovine) -- Novartis  |
| IGN-101 -- Igeneon   | insulin analogue -- Eli Lilly   |
| IK HIR02 -- Iketon   | Insulin Aspart -- Novo Nordisk  |
| IL-11 -- Genetics Institute/AHP                              | insulin detemir -- Novo Nordisk   |
| IL-13-PE38 -- NeoPharm                                       | insulin glargine -- Aventis   |
| IL-17 receptor -- Immunex                                    | insulin inhaled -- Inhale Therapeutics Systems, Alkermes  |
| IL-18BP -- Yeda  | insulin oral -- Inovax  |
| IL-1Hy1 -- Hyseq   | insulin, AeroDose -- AeroGen  |
| IL-1 $\beta$ -- Celltech                                     | insulin, AERx -- Aradigm  |
| IL-1 $\beta$ adjuvant -- Celltech                            | insulin, BEODAS -- Elan   |
| IL-2 -- Chiron   | insulin, Biphax -- Helix  |
| IL-2 + IL-12 -- Hoffman La-Roche                             | insulin, buccal -- Generec  |
| IL-6/sIL-6R fusion -- Hadasit                                | insulin, I2R -- Flemington  |
| IL-6R derivative -- Tosoh                                    | insulin, intranasal -- Bentley  |
| IL-7-Dap 389 fusion toxin -- Ligand                          | insulin, oral -- Nobex, Unigene   |
| IL-21 -- Novo Nordisk, ZymoGenetics                          | insulin, Orasome -- Endorec   |
| IM-862 -- Cytran   | insulin, ProMaxx -- Epic  |
| IMC-1C11 -- ImClone  | insulin, Quadrant -- Elan   |
| imiglucerase -- Genzyme                                      | insulin, recombinant -- Aventis   |
| Immune globulin intravenous (human) -- Hoffman La Roche      | insulin, Spiros -- Elan   |
| immune privilege factor -- Proneuron                         | insulin, Transfersome -- IDEA   |
| Immunocal -- Immunotec                                       | insulin, Zymo, recombinant -- Novo Nordisk  |
| Immunogene therapy -- Briana Bio-Tech                        | insulinotropin -- Scios   |
| Immunoliposomal 5-fluorodeoxyuridine-dipalmitate --          | Insulysin gene therapy --   |
| immunosuppressant vaccine -- Aixlie                          | integrin antagonists -- Merck   |
| immunotoxin -- Antisoma, NIH                                 | interferon (Alpha2) -- SRC VB VECTOR, Viragen, Dong-A, Hoffman La-Roche, Genentech  |
| ImmuRAIT-Re-188 -- Immunomedics                              | interferon -- BioMedicines, Human Genome Sciences   |
| imreg-1 -- Imreg   | interferon (Alfa-n3) -- Interferon Sciences Intl.   |
| infertility -- Johnson & Johnson, E-TRANS                    | interferon (Alpha), Biphax -- Helix   |
| Infliximab -- Centocor                                       |   |
| Influenza virus vaccine -- Aventis Pasteur, Protein Sciences |   |
| inhibin -- Biotech Australia, Human Therapeutics             |   |
| Inhibitory G protein gene therapy                            |   |
| INKP-2001 -- InKine  |   |
| Inolimomab -- Diaclone                                       |   |

FIG. 280

## 46/498

|   |   |
|---|---|
| interferon (Alpha)—Amgen, BioNative,        | IL-2/ diphtheria toxin -- Ligand              |
| Novartis, Genzyme Transgenics,              | Interleukin-3 -- Cangene                      |
| Hayashibara, Inhale Therapeutics            | Interleukin-4 -- Immunology Ventures,         |
| Systems, Medusa, Flamel, Dong-A,            | Sanofi Winthrop, Schering-Plough,             |
| GeneTrol, Nastech, Shantha,                 | Immunex/ Sanofi Winthrop, Bayer, Ono          |
| Wassermann, LG Chem, Sumitomo,              | interleukin-4 + TNF-Alpha -- NIH              |
| Aventis, Behring EGIS, Pepgen, Servier,     | interleukin-4 agonist -- Bayer                |
| Rhein Biotech,                              | interleukin-4 fusion toxin -- Ligand          |
| interferon (Alpha2A)                        | Interleukin-4 receptor -- Immunex, Immun      |
| interferon (Alpha2B) -- Enzon, Schering-    | Interleukin-6 -- Ajinomoto, Cangene, Yeda,    |
| Plough, Biogen, IDEA                        | Genetics Institute, Novartis                  |
| interferon (Alpha-N1) -- GlaxoSmithKline    | interleukin-6 fusion protein                  |
| interferon (beta) -- Rentschler, GeneTrol,  | interleukin-6 fusion toxin -- Ligand, Serono  |
| Meristem, Rhein Biotech, Toray, Yeda,       | interleukin-7 -- IC Innovations               |
| Daiichi, Mochida                            | interleukin-7 receptor -- Immunex             |
| interferon (Beta1A) -- Serono, Biogen       | interleukin-8 antagonists -- Kyowa            |
| interferon (beta1A),inhale -- Biogen        | Hakko/Millennium/Pfizer                       |
| interferon (β1b)-- Chiron                   | interleukin-9 antagonists -- Genaera          |
| interferon (tau)-- Pepgen                   | Interleukin-10 -- DNAX, Schering-Plough       |
| Interferon alfacon-1 -- Amgen               | Interleukin-10 gene therapy --                |
| Interferon alpha-2a vaccine                 | interleukin-12 -- Genetics Institute, Hoffman |
| Interferon Beta 1b -- Schering/Chiron,      | La-Roche                                      |
| InterMune                                   | interleukin-13 -- Sanofi                      |
| Interferon Gamma -- Boehringer Ingelheim,   | interleukin-13 antagonists -- AMRAD           |
| Sheffield, Rentschler, Hayashibara          | Interleukin-13-PE38QQR                        |
| interferon receptor , Type I -- Serono      | interleukin-15 -- Immunex                     |
| interferon(Gamma1B) -- Genentech            | interleukin-16 -- Research Corp               |
| Interferon-alpha-2b + ribavirin -- Biogen,  | interleukin-18 -- GlaxoSmithKline             |
| ICN   | Interleukin-18 binding protein -- Serono      |
| Interferon-alpha-2b gene therapy --         | Ior-P3 -- Center of Molecular Immunology      |
| Schering-Plough                             | IP-10 -- NIH                                  |
| Interferon-con1 gene therapy --             | IPF -- Metabolex                              |
| interleukin-1 antagonists -- Dompe          | IR-501 -- Immune Response                     |
| Interleukin-1 receptor antagonist -- Abbott | ISIS 9125 -- Isis Pharmaceuticals             |
| Bioresearch, Pharmacia                      | ISURF No. 1554 -- Millennium                  |
| Interleukin-1 receptor type I -- Immunex    | ISURF No. 1866 -- Iowa State Univer.          |
| interleukin-1 receptor Type II -- Immunex   | ITF-1697 -- Italfarmaco                       |
| Interleukin-1 trap -- Regeneron             | IxC 162 -- Ixion                              |
| Interleukin-1-alpha -- Immunex/Roche        | J 695 -- Cambridge Antibody Tech.,            |
| interleukin-2 -- SRC VB VECTOR,             | Genetics Inst., Knoll                         |
| Ajinomoto, Biomira, Chiron                  | Jagged + FGF -- Repair                        |

FIG. 28P

## 47/498

|  |   |
|--|---|
| JKC-362 -- Phoenix Pharmaceuticals             | leptin, 2nd-generation -- Amgen   |
| JTP-2942 -- Japan Tobacco                      | leridistim -- Pharmacia   |
| Juman monoclonal antibodies -- Medarex/Raven   | leuprolide, ProMaxx -- Epic   |
| K02 -- Axys Pharmaceuticals                    | leuprorelin, oral -- Unigene  |
| Keliximab -- IDEC                              | LeuTech -- Papatin  |
| Keyhole limpet haemocyanin                     | LEX 032 -- SuperGen   |
| KGF -- Amgen                                   | LiDEPT -- Novartis  |
| KM 871 -- Kyowa                                | Lintuzumab (anti-CD33 MAb) -- Protein Design Labs                       |
| KPI 135 -- Scios                               | lipase -- Altus Biologics   |
| KPI-022 -- Scios                               | lipid A vaccine -- EntreMed   |
| Kringle 5                                      | lipid-linked anchor Tech. -- ICRT, ID Biomedical                        |
| KSB 304  | liposome-CD4 Tech. -- Sheffield   |
| KSB-201 -- KS Biomedix                         | Listeria monocytogenes vaccine  |
| L 696418 -- Merck                              | LMB 1   |
| L 703801 -- Merck                              | LMB 7   |
| L1 -- Acorda                                   | LMB 9 -- Battelle Memorial Institute, NIH                               |
| L-761191 -- Merck                              | LM-CD45 -- Cantab Pharmaceuticals                                       |
| lactoferrin -- Meristem, Pharming, Agennix     | lovastatin -- Merck   |
| lactoferrin cardio -- Pharming                 | LSA-3   |
| LAG-3 -- Serono                                | LT- $\beta$ receptor -- Biogen  |
| LAIT -- GEMMA                                  | lung cancer vaccine -- Corixa   |
| LAK cell cytotoxin -- Arizona                  | lusupultide -- Scios  |
| lamellarins -- PharmaMar/University of Malaga  | L-Vax -- AVAX   |
| laminin A peptides -- NIH                      | LY 355455 -- Eli Lilly  |
| lanotepase -- Genetics Institute               | LY 366405 -- Eli Lilly  |
| laronidase -- BioMarin                         | LY-355101 -- Eli Lilly  |
| Lassa fever vaccine                            | Lyme disease DNA vaccine -- Vical/Aventis Pasteur                       |
| LCAT -- NIH                                    | Lyme disease vaccine -- Aquila  |
| LDP 01 -- Millennium                           | Biopharmaceuticals, Aventis, Pasteur, Symbicom, GlaxoSmithKline, Hyland |
| LDP 02 -- Millennium                           | Immuno, MedImmune   |
| Lecithinized superoxide dismutase -- Seikagaku | Lymphocytic choriomeningitis virus vaccine                              |
| LeIF adjuvant -- Corixa                        | lymphoma vaccine -- Biomira, Genitope                                   |
| leishmaniasis vaccine -- Corixa                | LYP18   |
| lenercept -- Hoffman La-Roche                  | lys plasminogen, recombinant  |
| Lenograstim -- Aventis, Chugai                 | Lysosomal storage disease gene therapy -- Avigen                        |
| lepirudin -- Aventis                           | lysostaphin -- Nutrition 21   |
| leptin -- Amgen, IC Innovations                |   |
| Leptin gene therapy -- Chiron Corporation      |   |

FIG. 28Q



48/498

|  |   |
|--|---|
| M 23 -- Gruenenthal                    | MEDI 507 -- BioTransplant                 |
| M1 monoclonal antibodies -- Acorda     | melanin concentrating hormone --          |
| Therapeutics                           | Neurocrine Biosciences                    |
| MA 16N7C2 -- Corvas Intl.              | melanocortins -- OMRF                     |
| malaria vaccine -- GlaxoSmithKline,    | Melanoma monoclonal antibodies -- Viragen |
| AdProTech, Antigenics, Apovia, Aventis | melanoma vaccine -- GlaxoSmithKline,      |
| Pasteur, Axis Genetics, Behringwerke,  | Akzo Nobel, Avant, Aventis Pasteur,       |
| CDCP, Chiron Vaccines, Genzyme         | Bavarian Nordic, Biovector, CancerVax,    |
| Transgenics, Hawaii, MedImmune, NIH,   | Genzyme Molecular Oncology, Humbolt,      |
| NYU, Oxxon, Roche/Saramane, Biotech    | ImClone Systems, Memorial, NYU, Oxxon     |
| Australia, Rx Tech                     | Melanoma vaccine Magevac -- Therion       |
| Malaria vaccine CDC/NIIMALVAC-1        | memory enhancers -- Scios                 |
| malaria vaccine, multicomponent        | meningococcal B vaccine -- Chiron         |
| mammaglobin -- Corixa                  | meningococcal vaccine -- CAMR             |
| mammastatin -- Biotherapeutics         | Meningococcal vaccine group B conjugate - |
| mannan-binding lectin -- NatlImmu      | - North American Vaccine                  |
| mannan-MUC1 -- Psiron                  | Meningococcal vaccine group B             |
| MAP 30                                 | recombinant -- BioChem Vaccines,          |
| Marinovir -- Phytera                   | Microscience                              |
| MARstem -- Maret                       | Meningococcal vaccine group Y conjugate - |
| MB-015 -- Mochida                      | - North American Vaccine                  |
| MBP -- ImmuLogic                       | Meningococcal vaccine groups A B and C    |
| MCI-028 -- Mitsubishi-Tokyo            | conjugate -- North American Vaccine       |
| MCIF -- Human Genome Sciences          | Mepolizumab -- GlaxoSmithKline            |
| MDC -- Advanced BioScience -- Akzo     | Metastatin -- EntreMed, Takeda            |
| Nobel, ICOS                            | Met-CkB7 -- Human Genome Sciences         |
| MDX 11 -- Medarex                      | met-enkephalin -- TNI                     |
| MDX 210 -- Medarex                     | METH-1 -- Human Genome Sciences           |
| MDX 22 -- Medarex                      | methioninase -- AntiCancer                |
| MDX 22                                 | Methionine lyase gene therapy --          |
| MDX 240 -- Medarex                     | AntiCancer                                |
| MDX 33                                 | Met-RANTES -- Genexa Biomedical,          |
| MDX 44 -- Medarex                      | Serono                                    |
| MDX 447 -- Medarex                     | Metreleptin                               |
| MDX H210 -- Medarex                    | Microtubule inhibitor MAb                 |
| MDX RA -- Houston BioTech., Medarex    | Immunogen/Abgenix                         |
| ME-104 -- Pharmexa                     | MGDF -- Kirin                             |
| Measles vaccine                        | MGV -- Progenics                          |
| Mecasermin -- Cephalon/Chiron, Chiron  | micrin -- Endocrine                       |
| MEDI 488 -- MedImmune                  | microplasmin -- ThromboGenics             |
| MEDI 500                               | MIF -- Genetics Institute                 |

FIG. 28R

## 49/498

|   |   |
|---|---|
| migration inhibitory factor -- NIH  | MAB 45-2D9- -- haematoporphyrin conjugate                       |
| Mim CD4.1 -- Xyte Therapies   | MAB 4B4   |
| mirostipen -- Human Genome Sciences   | MAB 4E3-CPA conjugate -- BCM Oncologia                          |
| Mitumomab (BEC-2) -- ImClone Systems, Merck KGaA  | MAB 4E3-daunorubicin conjugate                                  |
| MK 852 -- Merck   | MAB 50-6  |
| MLN 1202 (Anti-CCR2 monoclonal antibody) -- Millenium Pharmaceuticals   | MAB 50-61A -- Institut Pasteur                                  |
| Mobenakin -- NIS  | MAB 5A8 -- Biogen   |
| molgramostim -- Genetics Institute, Novartis  | MAB 791T/36-methotrexate conjugate                              |
| monoclonal antibodies -- Abgenix/Celltech, Immusol/ Medarex, Viragen/ Roslin Institute, Cambridge Antibody Tech./Elan | MAB 7c11.e8   |
| MAB 108 --  | MAB 7E11 C5-selenocystamine conjugate                           |
| MAB 10D5 --   | MAB 93KA9 -- Novartis   |
| MAB 14.18-interleukin-2 immunocytokine -- Lexigen   | MAB A5B7-cisplatin conjugate -- Biodynamics Research, Pharmacia |
| MAB 14G2a --  | MAB A5B7-I-131  |
| MAB 15A10 --  | MAB A7  |
| MAB 170 -- Biomira  | MAB A717 -- Exocell   |
| MAB 177Lu CC49 --   | MAB A7-zinostatin conjugate                                     |
| MAB 17F9  | MAB ABX-RB2 -- Abgenix  |
| MAB 1D7   | MAB ACA 11  |
| MAB 1F7 -- Immune Network   | MAB AFP-I-131 -- Immunomedics                                   |
| MAB 1H10-doxorubicin conjugate  | MAB AP1   |
| MAB 26-2F   | MAB AZ1   |
| MAB 2A11  | MAB B3-LysPE40 conjugate  |
| MAB 2E1 -- RW Johnson   | MAB B4 -- United Biomedical                                     |
| MAB 2F5   | MAB B43 Genistein-conjugate                                     |
| MAB 31.1 -- International BioImmune Systems   | MAB B43.13-Tc-99m -- Biomira                                    |
| MAB 32 -- Cambridge Antibody Tech., Peptech   | MAB B43-PAP conjugate   |
| MAB 323A3 -- Centocor   | MAB B4G7-gelonin conjugate                                      |
| MAB 3C5   | MAB BCM 43-daunorubicin conjugate -- BCM Oncologia              |
| MAB 3F12  | MAB BIS-1   |
| MAB 3F8   | MAB BMS 181170 -- Bristol-Myers Squibb                          |
| MAB 42/6  | MAB BR55-2  |
| MAB 425 -- Merck KGaA   | MAB BW494   |
| MAB 447-52D -- Merck Sharp & Dohme  | MAB C 242-DM1 conjugate -- ImmunoGen                            |
|   | MAB C242-PE conjugate   |
|   | MAB c30-6   |
|   | MAB CA208-cytorhodin-S conjugate -- Hoechst Japan               |
|   | MAB CC49 -- Enzon   |

FIG. 28S

## 50/498

|   |   |
|---|---|
| MAB ch14.18 --                          | MAB LL2-I-131 -- Immunomedics           |
| MAB CH14.18-GM-CSF fusion protein --    | MAB LL2-Y-90                            |
| Lexigen                                 | MAB LS2D617 -- Hybritech                |
| MAB chCE7                               | MAB LYM-1-gelonin conjugate             |
| MAB CI-137 -- AMRAD                     | MAB LYM-1-I-131                         |
| MAB cisplatin conjugate                 | MAB LYM-1-Y-90                          |
| MAB CLB-CD19                            | MAB LYM-2 -- Peregrine                  |
| MAB CLB-CD19v                           | MAB M195                                |
| MAB CLL-1 -- Peregrine                  | MAB M195-bismuth 213 conjugate --       |
| MAB CLL-1-GM-CSF conjugate              | Protein Design Labs                     |
| MAB CLL-1-IL-2 conjugate -- Peregrine   | MAB M195-gelonin conjugate              |
| MAB CLN IgG -- doxorubicin conjugates   | MAB M195-I-131                          |
| MAB conjugates -- Tanox                 | MAB M195-Y-90                           |
| MAB D612                                | MAB MA 33H1 -- Sanofi                   |
| MAB Dal B02                             | MAB MAD11                               |
| MAB DC101 -- ImClone                    | MAB MGB2                                |
| MAB EA 1 --                             | MAB MINT5                               |
| MAB EC708 -- Biovation                  | MAB MK2-23                              |
| MAB EP-5C7 -- Protein Design Labs       | MAB MOC31 ETA(252-613) conjugate        |
| MAB ERIC-1 -- ICRT                      | MAB MOC-31-In-111                       |
| MAB F105 gene therapy                   | MAB MOC-31-PE conjugate                 |
| MAB FC 2.15                             | MAB MR6 --                              |
| MAB G250 -- Centocor                    | MAB MRK-16 -- Aventis Pasteur           |
| MAB GA6                                 | MAB MS11G6                              |
| MAB GA733                               | MAB MX-DTPA BrE-3                       |
| MAB Gliomab-H -- Viventia Biotech       | MAB MY9                                 |
| MAB HB2-saporin conjugate               | MAB Nd2 -- Tosoh                        |
| MAB HD 37 --                            | MAB NG-1 -- Hygeia                      |
| MAB HD37-ricin chain-A conjugate        | MAB NM01 -- Nissin Food                 |
| MAB HNK20 -- Acambis                    | MAB OC 125                              |
| MAB huN901-DM1 conjugate --             | MAB OC 125-CMA conjugate                |
| ImmunoGen                               | MAB OKI-1 -- Ortho-McNeil               |
| MAB I-131 CC49 -- Corixa                | MAB OX52 -- Bioproducts for Science     |
| MAB ICO25                               | MAB PMA5                                |
| MAB ICR12-CPG2 conjugate                | MAB PR1                                 |
| MAB ICR-62                              | MAB prost 30                            |
| MAB IRac-ricin A conjugate              | MAB R-24                                |
| MAB K1                                  | MAB R-24 $\alpha$ Human GD3 -- Celltech |
| MAB KS1-4-methotrexate conjugate        | MAB RFB4-ricin chain A conjugate        |
| MAB L6 -- Bristol-Myers Squibb, Oncogen | MAB RFT5-ricin chain A conjugate        |
| MAB LiCO 16-88                          | MAB SC 1                                |

FIG. 28T

## 51/498

|  |  |
|--|--|
| MAb SM-3 -- ICRT                       | Muc-1 vaccine -- Corixa                    |
| MAb SMART 1D10 -- Protein Design Labs  | mucosal tolerance -- Aberdeen              |
| MAb SMART ABL 364 -- Novartis          | mullerian inhibiting subst                 |
| MAb SN6f                               | muplestim -- Genetics Institute, Novartis, |
| MAb SN6f-deglycosylated ricin A chain  | DSM Anti-Infectives                        |
| conjugate --                           | murine MAb -- KS Biomedix                  |
| MAb SN6j                               | Mutant somatropin -- JCR Pharmaceutical    |
| MAb SN7-ricin chain A conjugate        | MV 833 -- Toagosei                         |
| MAb T101-Y-90 conjugate -- Hybritech   | Mycoplasma pulmonis vaccine                |
| MAb T-88 -- Chiron                     | Mycoprex -- XOMA                           |
| MAb TB94 -- Cancer ImmunoBiology       | myeloperoxidase -- Henogen                 |
| MAb TEC 11                             | myostatin -- Genetics Institute            |
| MAb TES-23 -- Chugai                   | Nacolomab tafenatox -- Pharmacia           |
| MAb TM31 -- Avant                      | Nagrecor -- Scios                          |
| MAb TNT-1 -- Cambridge Antibody Tech., | nagrestipen -- British Biotech             |
| Peregrine                              | NAP-5 -- Corvas Intl.                      |
| MAb TNT-3                              | NAPc2 -- Corvas Intl.                      |
| MAb TNT-3 -- IL2 fusion protein --     | nartograstim -- Kyowa                      |
| MAb TP3-A $\alpha$ -211                | Natalizumab -- Protein Design Labs         |
| MAb TP3-PAP conjugate --               | Nateplase -- NIH, Nihon Schering           |
| MAb UJ13A -- ICRT                      | nateplase -- Schering AG                   |
| MAb UN3                                | NBI-3001 -- Neurocrine Biosci.             |
| MAb ZME-018-gelonin conjugate          | NBI-5788 -- Neurocrine Biosci.             |
| MAb-BC2 -- GlaxoSmithKline             | NBI-6024 -- Neurocrine Biosci.             |
| MAb-DM1 conjugate -- ImmunoGen         | Nef inhibitors -- BRI                      |
| MAb-ricin-chain-A conjugate -- XOMA    | Neisseria gonorrhoea vaccine -- Antex      |
| MAb-temoporfin conjugates              | Biologics                                  |
| Monopharm C -- Viventia Biotech        | Neomycin B-arginine conjugate              |
| monteplase -- Eisai                    | Nerelimomab -- Chiron                      |
| montirelin hydrate -- Gruenenthal      | Nerve growth factor -- Amgen -- Chiron,    |
| moroctocog alfa -- Genetics Institute  | Genentech                                  |
| Moroctocog-alfa -- Pharmacia           | Nerve growth factor gene therapy           |
| MP 4                                   | nesiritide citrate -- Scios                |
| MP-121 -- Biopharm                     | neuregulin-2 -- CeNeS                      |
| MP-52 -- Biopharm                      | neurocan -- NYU                            |
| MRA -- Chugai                          | neuronal delivery system -- CAMR           |
| MS 28168 -- Mitsui Chemicals, Nihon    | Neurophil inhibitory Factor -- Corvas      |
| Schering                               | Neuroprotective vaccine -- University of   |
| MSH fusion toxin -- Ligand             | Auckland                                   |
| MSI-99 -- Genaera                      | neurotrophic chimaeras -- Regeneron        |
| MT 201 -- Micromet                     | neurotrophic factor -- NsGene, CereMedix   |

FIG. 28U

## 52/498

|                                    |  |
|------------------------------------|--|
| NeuroVax -- Immune Response        | Oncophage -- Antigenics                    |
| neurturin -- Genentech             | Oncostatin M -- Bristol-Myers Squibb       |
| neutral endopeptidase -- Genentech | OncoVax-CL -- Jenner Biotherapies          |
| NGF enhancers -- NeuroSearch       | OncoVax-P -- Jenner Biotherapies           |
| NHL vaccine -- Large Scale Biology | onercept -- Yeda                           |
| NIP45 -- Boston Life Sciences      | onychomycosis vaccine -- Boehringer        |
| NKI-B20                            | Ingelheim                                  |
| NM 01 -- Nissin Food               | opebecan -- XOMA                           |
| NMI-139 -- NitroMed                | opioids -- Arizona                         |
| NMMP -- Genetics Institute         | Oprelvekin -- Genetics Institute           |
| NN-2211 -- Novo Nordisk            | Oregovomab -- AltaRex                      |
| Noggin -- Regeneron                | Org-33408 b-- Akzo Nobel                   |
| Nonacog alfa                       | Orolip DP -- EpiCept                       |
| Norelin -- Biostar                 | oryzacystatin                              |
| Norwalk virus vaccine              | OSA peptides -- GenSci Regeneration        |
| NRLU 10 -- NeoRx                   | osteoblast-cadherin GF -- Pharis           |
| NRLU 10 PE -- NeoRx                | Osteocalcin-thymidine kinase gene therapy  |
| NT-3 -- Regeneron                  | osteogenic protein -- Curis                |
| NT-4/5 -- Genentech                | osteopontin -- OraPharma                   |
| NU 3056                            | osteoporosis peptides -- Integra, Telios   |
| NU 3076                            | osteoprotegerin -- Amgen, SnowBrand        |
| NX 1838 -- Gilead Sciences         | otitis media vaccines -- Antex Biologics   |
| NY ESO-1/CAG-3 antigen -- NIH      | ovarian cancer -- University of Alabama    |
| NYVAC-7 -- Aventis Pasteur         | OX40-IgG fusion protein -- Cantab, Xenova  |
| NZ-1002 -- Novazyme                | P 246 -- Diatide                           |
| obesity therapy -- Nobex           | P 30 -- Alfacell                           |
| OC 10426 -- Ontogen                | p1025 -- Active Biotech                    |
| OC 144093 -- Ontogen               | P-113 <sup>^</sup> -- Demegen              |
| OCIF -- Sankyo                     | P-16 peptide -- Transition Therapeutics    |
| Oct-43 -- Otsuka                   | p43 -- Ramot                               |
| Odulimomab -- Immunotech           | P-50 peptide -- Transition Therapeutics    |
| OK PSA - liposomal                 | p53 + RAS vaccine -- NIH, NCI              |
| OKT3-gamma-1-ala-ala               | PACAP(1-27) analogue                       |
| OM 991                             | paediatric vaccines -- Chiron              |
| OM 992                             | Pafase -- ICOS                             |
| Omalizumab -- Genentech            | PAGE-4 plasmid DNA -- IDEC                 |
| oncoimmunin-L -- NIH               | PAI-2 -- Biotech Australia, Human          |
| Oncolysin B -- ImmunoGen           | Therapeutics                               |
| Oncolysin CD6 -- ImmunoGen         | Palifermin (keratinocyte growth factor) -- |
| Oncolysin M -- ImmunoGen           | Amgen                                      |
| Oncolysin S -- ImmunoGen           | Palivizumab -- MedImmune                   |

FIG. 28V

## 53/498

|   |   |
|---|---|
| PAM 4 -- Merck                            | PEG-uricase -- Mountain View              |
| pamiteplase -- Yamanouchi                 | Pegvisomant -- Genentech                  |
| pancreatin, Minitabs -- Eurand            | PEGylated proteins, PolyMASC -- Valentis  |
| Pangen -- Fournier                        | PEGylated recombinant native human leptin |
| Pantarin -- Selective Genetics            | -- Roche                                  |
| Parainfluenza virus vaccine -- Pharmacia, | Pemtumomab                                |
| Pierre Fabre                              | Penetratin -- Cyclacel                    |
| paraoxanase -- Esperion                   | Pepscan -- Antisoma                       |
| parathyroid hormone -- Abiogen, Korea     | peptide G -- Peptech, ICRT                |
| Green Cross                               | peptide vaccine -- NIH, NCI               |
| Parathyroid hormone (1-34) --             | Pexelizumab                               |
| Chugai/Suntory                            | pexiganan acetate -- Genaera              |
| Parkinson's disease gene therapy -- Cell  | Pharmaprojects No. 3179 -- NYU            |
| Genesys/ Ceregene                         | Pharmaprojects No. 3390 -- Ernest Orlando |
| Parvovirus vaccine -- MedImmune           | Pharmaprojects No. 3417 -- Sumitomo       |
| PCP-Scan -- Immunomedics                  | Pharmaprojects No. 3777 -- Acambis        |
| PDGF -- Chiron                            | Pharmaprojects No. 4209 -- XOMA           |
| PDGF cocktail -- Theratechnologies        | Pharmaprojects No. 4349 -- Baxter Intl.   |
| peanut allergy therapy -- Dynavax         | Pharmaprojects No. 4651                   |
| PEG anti-ICAM MAb -- Boehringer           | Pharmaprojects No. 4915 -- Avanir         |
| Ingelheim                                 | Pharmaprojects No. 5156 -- Rhizogenics    |
| PEG asparaginase -- Enzon                 | Pharmaprojects No. 5200 -- Pfizer         |
| PEG glucocerebrosidase                    | Pharmaprojects No. 5215 -- Origene        |
| PEG hirudin -- Knoll                      | Pharmaprojects No. 5216 -- Origene        |
| PEG interferon-alpha-2a -- Roche          | Pharmaprojects No. 5218 -- Origene        |
| PEG interferon-alpha-2b + ribavirin --    | Pharmaprojects No. 5267 -- ML             |
| Biogen, Enzon, ICN Pharmaceuticals,       | Laboratories                              |
| Schering-Plough                           | Pharmaprojects No. 5373 -- MorphoSys      |
| PEG MAb A5B7 --                           | Pharmaprojects No. 5493 -- Metabolex      |
| Pegacaristim -- Amgen -- Kirin Brewery -- | Pharmaprojects No. 5707 -- Genentech      |
| ZymoGenetics                              | Pharmaprojects No. 5728 -- Autogen        |
| Pegaldesleukin -- Research Corp           | Pharmaprojects No. 5733 -- BioMarin       |
| pegaspargase -- Enzon                     | Pharmaprojects No. 5757 -- NIH            |
| pegfilgrastim -- Amgen                    | Pharmaprojects No. 5765 -- Gryphon        |
| PEG-interferon Alpha -- Viragen           | Pharmaprojects No. 5830 -- AntiCancer     |
| PEG-interferon Alpha 2A -- Hoffman La-    | Pharmaprojects No. 5839 -- Dyax           |
| Roche                                     | Pharmaprojects No. 5849 -- Johnson &      |
| PEG-interferon Alpha 2B -- Schering-      | Johnson                                   |
| Plough                                    | Pharmaprojects No. 5860 -- Mitsubishi-    |
| PEG-r-hirudin -- Abbott                   | Tokyo                                     |
| PEG-rHuMGDF -- Amgen                      |   |

FIG. 28W

## 54/498

|   |   |
|---|---|
| Pharmaprojects No. 5869 -- Oxford GlycoSciences   | Plasminogen activators -- Abbott Laboratories, American Home Products, Boehringer Mannheim, Chiron  |
| Pharmaprojects No. 5883 -- Asahi Brewery  | Corporation, DuPont Pharmaceuticals, Eli Lilly, Shionogi, Genentech, Genetics Institute, GlaxoSmithKline, Hemispherx Biopharma, Merck & Co, Novartis, Pharmacia Corporation, Wakamoto, Yeda |
| Pharmaprojects No. 5947 -- StressGen  | plasminogen-related peptides -- Bio-Tech. General/MGH   |
| Pharmaprojects No. 5961 -- Theratechnologies  | platelet factor 4 -- RepliGen   |
| Pharmaprojects No. 5962 -- NIH  | Platelet-derived growth factor -- Amgen -- ZymoGenetics   |
| Pharmaprojects No. 5966 -- NIH  | plusonemin-- Hayashibara  |
| Pharmaprojects No. 5994 -- Pharming   | PMD-2850 -- Protherics  |
| Pharmaprojects No. 5995 -- Pharming   | Pneumococcal vaccine -- Antex Biologics, Aventis Pasteur  |
| Pharmaprojects No. 6023 -- IMMUCON  | Pneumococcal vaccine intranasal -- BioChem Vaccines/Biovector   |
| Pharmaprojects No. 6063 -- Cytoclonal   | PR1A3   |
| Pharmaprojects No. 6073 -- SIDDCO   | PR-39   |
| Pharmaprojects No. 6115 -- Genzyme  | pralmorelin -- Kaken  |
| Pharmaprojects No. 6227 -- NIH  | Pretarget-Lymphoma -- NeoRx   |
| Pharmaprojects No. 6230 -- NIH  | Priliximab -- Centocor  |
| Pharmaprojects No. 6236 -- NIH  | PRO 140 -- Progenics  |
| Pharmaprojects No. 6243 -- NIH  | PRO 2000 -- Procept   |
| Pharmaprojects No. 6244 -- NIH  | PRO 367 -- Progenics  |
| Pharmaprojects No. 6281 -- Senetek  | PRO 542 -- Progenics  |
| Pharmaprojects No. 6365 -- NIH  | pro-Apo A-I -- Esperion   |
| Pharmaprojects No. 6368 -- NIH  | prolactin -- Genzyme  |
| Pharmaprojects No. 6373 -- NIH  | Prosaptide TX14(A) -- Bio-Tech. General   |
| Pharmaprojects No. 6408 -- Pan Pacific  | prostate cancer antibodies -- Immunex, UroCor   |
| Pharmaprojects No. 6410 -- Athersys   | prostate cancer antibody therapy -- Genentech/UroGenesys, Genotherapeutics  |
| Pharmaprojects No. 6421 -- Oxford GlycoSciences   | prostate cancer immunotherapeutics -- The PSMA Development Company  |
| Pharmaprojects No. 6522 -- Maxygen  | prostate cancer vaccine -- Aventis Pasteur, Zonagen, Corixa, Dendreon, Jenner   |
| Pharmaprojects No. 6523 -- Pharis   | Biotherapies, Therion Biologics   |
| Pharmaprojects No. 6538 -- Maxygen  |   |
| Pharmaprojects No. 6554 -- APALEXO  |   |
| Pharmaprojects No. 6560 -- Ardana   |   |
| Pharmaprojects No. 6562 -- Bayer  |   |
| Pharmaprojects No. 6569 -- Eos  |   |
| Phenoxazine   |   |
| Phenylase -- Ibex   |   |
| Pigment epithelium derived factor -- plasminogen activator inhibitor-1, recombinant -- DuPont Pharmaceuticals |   |

FIG. 28X

## 55/498

|  |  |
|--|--|
| prostate-specific antigen -- EntreMed        | RD 62198                                     |
| protein A -- RepliGen                        | rDnase -- Genentech                          |
| protein adhesives -- Enzon                   | RDP-58 -- SangStat                           |
| protein C -- Baxter Intl., PPL Therapeutics, | RecepTox-Fce -- Keryx                        |
| ZymoGenetics                                 | RecepTox-GnRH -- Keryx, MTR                  |
| protein C activator -- Gilead Sciences       | Technologies                                 |
| protein kinase R antags -- NIH               | RecepTox-MBP -- Keryx, MTR                   |
| protirelin -- Takeda                         | Technologies                                 |
| protocadherin 2 -- Caprion                   | recFSH -- Akzo Nobel, Organon                |
| Pro-urokinase -- Abbott, Bristol-Myers       | REGA 3G12                                    |
| Squibb, Dainippon, Tosoh -- Welfide          | Regavirumab -- Teijin                        |
| P-selectin glycoprotein ligand-1 -- Genetics | relaxin -- Connetics Corp                    |
| Institute                                    | Renal cancer vaccine -- Macropharm           |
| pseudomonal infections -- InterMune          | repifermin -- Human Genome Sciences          |
| Pseudomonas vaccine -- Cytovax               | Respiratory syncytial virus PFP-2 vaccine -- |
| PSGL-Ig -- American Home Products            | Wyeth-Lederle                                |
| PSP-94 -- Procyon                            | Respiratory syncytial virus vaccine --       |
| PTH 1-34 -- Nobex                            | GlaxoSmithKline, Pharmacia, Pierre Fabre     |
| Quilimmune-M -- Antigenics                   | Respiratory syncytial virus vaccine          |
| R 744 -- Roche                               | inactivated                                  |
| R 101933                                     | Respiratory syncytial virus-parainfluenza    |
| R 125224 -- Sankyo                           | virus vaccine -- Aventis Pasteur,            |
| RA therapy -- Cardion                        | Pharmacia                                    |
| Rabies vaccine recombinant -- Aventis        | Retepase -- Boehringer Mannheim,             |
| Pasteur, BioChem Vaccines, Kaketsuken        | Hoffman La-Roche                             |
| Pharmaceuticals                              | Retropep -- Retroscreen                      |
| RadioTheraCIM -- YM BioSciences              | RFB4 (dsFv) PE38                             |
| Ramot project No. 1315 -- Ramot              | RFI 641 -- American Home Products            |
| Ramot project No. K-734A -- Ramot            | RFTS -- UAB Research Foundation              |
| Ramot project No. K-734B -- Ramot            | RG 12986 -- Aventis Pasteur                  |
| Ranibizumab (Anti-VEGF fragment) --          | RG 83852 -- Aventis Pasteur                  |
| Genentech                                    | RG-1059 -- RepliGen                          |
| RANK -- Immunex                              | rGCR -- NIH                                  |
| ranpirnase -- Alfacell                       | rGLP-1 -- Restoragen                         |
| ranpirnase-anti-CD22 MAb -- Alfacell         | rGRF -- Restoragen                           |
| RANTES inhibitor -- Milan                    | rh Insulin -- Eli Lilly                      |
| RAPID drug delivery systems -- ARIAD         | RHAMM targeting peptides -- Cangene          |
| rasburicase -- Sanofi                        | rHb1.1 -- Baxter Intl.                       |
| rBPI-21, topical -- XOMA                     | rhCC10 -- Claragen                           |
| RC 529 -- Corixa                             | rhCG -- Serono                               |
| rCFTR -- Genzyme Transgenics                 | Rheumatoid arthritis gene therapy            |

FIG. 28Y



## 56/498

|   |  |
|---|--|
| Rheumatoid arthritis vaccine -- Veterans Affairs Medical Center | SB RA 31012 --   |
| rhLH -- Serono  | SC 56929 -- Pharmacia                                    |
| Ribozyme gene therapy -- Genset                                 | SCA binding proteins -- Curis, Enzon                     |
| Rickettsial vaccine recombinant                                 | scFv(14E1)-ETA Berlex Laboratories, Schering AG          |
| RIGScan CR -- Neoprobe  | ScFv(FRP5)-ETA --  |
| RIP-3 -- Rigel  | ScFv6C6-PE40 --  |
| Rituximab -- Genentech  | SCH 55700 -- Celltech                                    |
| RK-0202 -- RxKinetix  | Schistosomiasis vaccine -- Glaxo Wellcome/Medeva, Brazil |
| RLT peptide -- Esperion   | SCPF -- Advanced Tissue Sciences                         |
| rM/NEI -- IVAX  | scuPA-suPAR complex -- Hadasit                           |
| rmCRP -- Immtech  | SD-9427 -- Pharmacia                                     |
| RN-1001 -- Renovo   | SDF-1 -- Ono   |
| RN-3 -- Renovo  | SDZ 215918 -- Novartis                                   |
| RNAse conjugate -- Immunomedics                                 | SDZ 280125 -- Novartis                                   |
| RO 631908 -- Roche  | SDZ 89104 -- Novartis                                    |
| Rotavirus vaccine -- Merck                                      | SDZ ABL 364 -- Novartis                                  |
| RP 431 -- DuPont Pharmaceuticals                                | SDZ MMA 383 -- Novartis                                  |
| RP-128 -- Resolution  | Secretin -- Ferring, Repligen                            |
| RPE65 gene therapy --   | serine protease inhbs -- Pharis                          |
| RPR 110173 -- Aventis Pasteur                                   | sermorelin acetate -- Serono                             |
| RPR 115135 -- Aventis Pasteur                                   | SERP-1 -- Viron  |
| RPR 116258A -- Aventis Pasteur                                  | sertenef -- Dainippon                                    |
| rPSGL-Ig -- American Home Products                              | serum albumin, Recombinant human -- Aventis Behring      |
| r-SPC surfactant -- Byk Gulden                                  | serum-derived factor -- Hadasit                          |
| RSV antibody -- Medimmune                                       | Sevirumab -- Novartis                                    |
| Ruplizumab -- Biogen  | SGN 14 -- Seattle Genetics                               |
| rV-HER-2/neu -- Therion Biologics                               | SGN 15 -- Seattle Genetics                               |
| SA 1042 -- Sankyo   | SGN 17/19 -- Seattle Genetics                            |
| sacrosidase -- Orphan Medical                                   | SGN 30 -- Seattle Genetics                               |
| Sant 7  | SGN-10 -- Seattle Genetics                               |
| Sargramostim -- Immunex   | SGN-11 -- Seattle Genetics                               |
| saruplase -- Gruenenthal  | SH 306 -- DuPont Pharmaceuticals                         |
| Satumomab -- Cytogen  | Shanvac-B -- Shantha                                     |
| SB 1 -- COR Therapeutics  | Shigella flexneri vaccine -- Avant, Acambis, Novavax     |
| SB 207448 -- GlaxoSmithKline                                    | Shigella sonnei vaccine --                               |
| SB 208651 -- GlaxoSmithKline                                    | sICAM-1 -- Boehringer Ingelheim                          |
| SB 240683 -- GlaxoSmithKline                                    | Silteplase -- Genzyme                                    |
| SB 249415 -- GlaxoSmithKline                                    |  |
| SB 249417 -- GlaxoSmithKline                                    |  |
| SB 6 -- COR Therapeutics  |  |

FIG. 28Z

## 57/498

|  |  |
|--|--|
| SIV vaccine -- Endocon, Institut Pasteur       | Staphylococcus aureus vaccine conjugate --   |
| SK 896 -- Sanwa Kagaku Kenkyusho               | Nabi   |
| SK-827 -- Sanwa Kagaku Kenkyusho               | Staphylococcus therapy -- Tripep             |
| Skeletex -- CellFactors                        | Staphylokinase -- Biovation, Prothera,       |
| SKF 106160 -- GlaxoSmithKline                  | Thrombogenetics                              |
| S-nitroso-AR545C --                            | Streptococcal A vaccine -- M6                |
| SNTP -- Active Biotech                         | Pharmaceuticals, North American Vaccine      |
| somatomedin-1 -- GroPep, Mitsubishi-           | Streptococcal B vaccine -- Microscience      |
| Tokyo, NIH                                     | Streptococcal B vaccine recombinant --       |
| somatomedin-1 carrier protein -- Insmed        | Biochem Vaccines                             |
| somatostatin -- Ferring                        | Streptococcus pyogenes vaccine               |
| Somatotropin/                                  | STRL-33 -- NIH                               |
| Human Growth Hormone -- Bio-Tech.              | Subalin -- SRC VB VECTOR                     |
| General, Eli Lilly                             | SUIS -- United Biomedical                    |
| somatropin -- Bio-Tech. General, Alkermes,     | SUIS-LHRH -- United Biomedical               |
| ProLease, Aventis Behring, Biovector,          | SUN-E3001 -- Suntory                         |
| Cangene, Dong-A, Eli Lilly, Emisphere,         | super high affinity monoclonal antibodies -- |
| Enact, Genentech, Genzyme Transgenics,         | YM BioSciences                               |
| Grandis/InfiMed, CSL, InfiMed, MacroMed,       | Superoxide dismutase -- Chiron, Enzon,       |
| Novartis, Novo Nordisk, Pharmacia              | Ube Industries, Bio-Tech, Yeda               |
| Serono, TranXenoGen                            | superoxide dismutase-2 -- OXIS               |
| somatropin derivative -- Schering AG           | suppressin -- UAB Research Foundation        |
| somatropin, AIR -- Eli Lilly                   | SY-161-P5 -- ThromboGenics                   |
| Somatropin, inhaled -- Eli Lilly/Alkermes      | SY-162 -- ThromboGenics                      |
| somatropin, Kabi -- Pharmacia                  | Systemic lupus erythematosus vaccine --      |
| somatropin, Orasome -- Novo Nordisk            | MedClone/VivoRx                              |
| Sonermin -- Dainippon Pharmaceutical           | T cell receptor peptides -- Xoma             |
| SP(V5.2)C -- Supertek                          | T cell receptor peptide vaccine              |
| SPf66  | T4N5 liposomes -- AGI Dermatics              |
| sphingomyelinase -- Genzyme                    | TACI, soluble -- ZymoGenetics                |
| SR 29001 -- Sanofi                             | targeted apoptosis -- Antisoma               |
| SR 41476 -- Sanofi                             | tasonermin -- Boehringer Ingelheim           |
| SR-29001 -- Sanofi                             | TASP   |
| SS1(dsFV)-PE38 -- NeoPharm                     | TASP-V                                       |
| $\beta$ 2 microglobulin -- Avidex              | Tat peptide analogues -- NIH                 |
| $\beta$ 2-microglobulin fusion proteins -- NIH | TBP I -- Yeda                                |
| $\beta$ -amyloid peptides -- CeNeS             | TBP II                                       |
| $\beta$ -defensin -- Pharis                    | TBV25H -- NIH                                |
| Staphylococcus aureus infections --            | Tc 99m ior cea1 -- Center of Molecular       |
| Inhibitex/ZLB                                  | Immunology                                   |
|  | Tc 99m P 748 -- Diatide                      |

FIG. 28AA

## 58/498

|  |   |
|--|---|
| Tc 99m votumumab -- Intracell                                      | Tissue factor -- Genentech  |
| Tc-99m rh-Annexin V -- Theseus Imaging                             | Tissue factor pathway inhibitor                                       |
| teceleukin -- Biogen   | TJN-135 -- Tsumura  |
| tenecteplase -- Genentech  | TM 27 -- Avant  |
| Teriparatide -- Armour Pharmaceuticals,<br>Asahi Kasei, Eli Lilly  | TM 29 -- Avant  |
| terlipressin -- Ferring  | TMC-151 -- Tanabe Seiyaku   |
| testisin -- AMRAD  | TNF tumour necrosis factor -- Asahi Kasei                             |
| Tetrafibricin -- Roche   | TNF Alpha -- CytImmune  |
| TFPI -- EntreMed   | TNF antibody -- Johnson & Johnson                                     |
| tgD-IL-2 -- Takeda   | TNF binding protein -- Amgen  |
| TGF-Alpha -- ZymoGenetics  | TNF degradation product -- Oncotech                                   |
| TGF- $\beta$ -- Kolon  | TNF receptor -- Immunex   |
| TGF- $\beta$ 2 -- Insmad   | TNF receptor 1, soluble -- Amgen                                      |
| TGF- $\beta$ 3 -- OSI  | TNF Tumour necrosis factor-alpha -- Asahi<br>Kasei, Genetech, Mochida |
| Thalassaemia gene therapy -- Crucell                               | TNF-Alpha inhibitor -- Tripep   |
| TheraCIM-h-R3 -- Center of Molecular<br>Immunology, YM BioSciences | TNFR:Fc gene therapy -- Targeted Genetics                             |
| Theradigm-HBV -- Epimmune  | TNF-SAM2  |
| Theradigm-HPV -- Epimmune  | Tolerimab -- Innogenetics   |
| Theradigm-malaria -- Epimmune                                      | Toxoplasma gondii vaccine --<br>GlaxoSmithKline                       |
| Theradigm-melanoma -- Epimmune                                     | TP 9201 -- Telios   |
| TheraFab -- Antisoma   | TP10 -- Avant   |
| ThGRF 1-29 -- Theratechnologies                                    | TP20 -- Avant   |
| ThGRF 1-44 -- Theratechnologies                                    | tPA -- Centocor   |
| Thrombin receptor activating peptide --<br>Abbott                  | trafermin -- Scios  |
| thrombomodulin -- Iowa, Novocastra                                 | TRAIL/Apo2L -- Immunex  |
| Thrombopoietin -- Dragon Pharmaceuticals,<br>Genentech             | TRAIL-R1 MAb -- Cambridge Antibody<br>Technologies                    |
| thrombopoietin, Pliva -- Recepton                                  | transferrin-binding proteins -- CAMR                                  |
| Thrombospondin 2 --  | Transforming growth factor-beta-1 --<br>Genentech                     |
| thrombostatin -- Thromgen  | transport protein -- Genesis  |
| thymalfasin -- SciClone  | Trastuzumab -- Genetech   |
| thymocartin -- Gedeon Richter                                      | TRH -- Ferring  |
| thymosin Alpha1 -- NIH   | Triabin -- Schering AG  |
| thyroid stimulating hormone -- Genzyme                             | Triconal  |
| tICAM-1 -- Bayer   | Triflavin   |
| Tick anticoagulant peptide -- Merck                                | troponin I -- Boston Life Sciences                                    |
| TIF -- Xoma  | TRP-2 <sup>A</sup> -- NIH   |
| Tifacogin -- Chiron, NIS, Pharmacia                                | trypsin inhibitor -- Mochida  |

FIG. 28BB

## 59/498

|   |  |
|---|--|
| TSP-1 gene therapy --                         | Vascular endothelial growth factors -- R&D |
| TT-232  | Systems                                    |
| TTS-CD2 -- Active Biotech                     | vascular targeting agents -- Peregrine     |
| Tuberculosis vaccine -- Aventis Pasteur,      | vasopermeation enhancement agents --       |
| Genesis                                       | Peregrine                                  |
| Tumor Targeted Superantigens -- Active        | vasostatin -- NIH                          |
| Biotech -- Pharmacia                          | VCL -- Bio-Tech. General                   |
| tumour vaccines -- PhotoCure                  | VEGF -- Genentech, Scios                   |
| tumour-activated prodrug antibody             | VEGF inhibitor -- Chugai                   |
| conjugates -- Millennium/ImmunoGen            | VEGF-2 -- Human Genome Sciences            |
| tumstatin -- ILEX                             | VEGF-Trap -- Regeneron                     |
| Tuvirumab -- Novartis                         | viscumin, recombinant -- Madaus            |
| TV-4710 -- Teva                               | Vitaxin                                    |
| TWEAK receptor -- Immunex                     | Vitrage -- ISTA Pharmaceuticals            |
| TXU-PAP                                       | West Nile virus vaccine -- Bavarian Nordic |
| TY-10721 -- TOA Eiyo                          | WP 652                                     |
| Type I diabetes vaccine -- Research Corp      | WT1 vaccine -- Corixa                      |
| Typhoid vaccine CVD 908                       | WX-293 -- Willex BioTech.                  |
| U 143677 -- Pharmacia                         | WX-360 -- Willex BioTech.                  |
| U 81749 -- Pharmacia                          | WX-UK1 -- Willex BioTech.                  |
| UA 1248 -- Arizona                            | XMP-500 -- XOMA                            |
| UGIF -- Sheffield                             | XomaZyme-791 -- XOMA                       |
| UIC 2   | XTL 001 -- XTL Biopharmaceuticals          |
| UK 101  | XTL 002 -- XTL Biopharmaceuticals          |
| UK-279276 -- Corvas Intl.                     | yeast delivery system -- GlobelImmune      |
| urodilatin -- Pharis                          | Yersinia pestis vaccine                    |
| urofollitrophin -- Serono                     | YIGSR-Stealth -- Johnson & Johnson         |
| Urokinase -- Abbott                           | Yisum Project No. D-0460 -- Yisum          |
| uteroferin-- Pepgen                           | YM 207 -- Yamanouchi                       |
| V 20 -- GLYCODESIGN                           | YM 337 -- Protein Design Labs              |
| V2 vasopressin receptor gene therapy          | Yttrium-90 labelled biotin                 |
| vaccines -- Active Biotech                    | Yttrium-90-labeled anti-CEA MAb T84.66 --  |
| Varicella zoster glycoprotein vaccine --      | ZD 0490 -- AstraZeneca                     |
| Research Corporation Technologies             | ziconotide -- Elan                         |
| Varicella zoster virus vaccine live -- Cantab | ZK 157138 -- Berlex Laboratories           |
| Pharmaceuticals                               | Zolimomab aritox                           |
| Vascular endothelial growth factor --         | Zorcell -- Immune Response                 |
| Genentech, University of California           | ZRXL peptides -- Novartis                  |

FIG. 28CC